

Specification

Proteinase K, recombinant

A7932

Solubility:	clear, colorless (20 mg/ml)
Product Code:	A7932
Product Name:	Proteinase K, recombinant
Short Description:	delivery form: lyophilized
Specifications:	DNases/RNases: not detectable Activity/mg: min. 30 mAnsonU Appearance: white lyophilizate Bioburden: ≤ 125 cfu/g DNA (threshold): ≤ 10 pg/mg enzyme
Hazard pictograms	 
WGK:	1
Storage:	2 - 8°C
Signal Word:	Danger
GHS Symbols:	GHS07 GHS08
H Phrases:	H315 H319 H334 H335
P Phrases:	P261 P284 P305+P351+P338 P342+P311 P405 P501

AppliChem GmbH

Ottoweg 4 • D-64291 Darmstadt • Phone +49 6151 9357 0 • Fax +49 6151 9357 11 • info.de@itwreagents.com • www.itwreagents.com
 CEO Joan Roget • Commerzbank Darmstadt • Bank 508 800 50 • Account 0186989900 IBAN DE24 5088 0050 0186 9899 00 • Swiftcode DRESDEFF508 • Finanzamt Darmstadt 0072250192 • Register court Darmstadt HRB Nr. 7340

Specification

Proteinase K, recombinant

A7932

M:	28.8 kD
CAS:	39450-01-6
EINECS:	254-457-8
CS:	35079090
Index Nr.:	647-014-00-9

from *Tritirachium album*

AppliChem GmbH

Ottoweg 4 • D-64291 Darmstadt • Phone +49 6151 9357 0 • Fax +49 6151 9357 11 • info.de@itwreagents.com • www.itwreagents.com
CEO Joan Roget • Commerzbank Darmstadt • Bank 508 800 50 • Account 0186989900 IBAN DE24 5088 0050 0186 9899 00 • Swiftcode DRESDEFF508 • Finanzamt Darmstadt 0072250192 • Register court Darmstadt HRB Nr. 7340