

Certificate of Conformity

VWR Catalogue Number	215-1908
Description	VWR® Rectangular Bottles, Wide Neck with Screw Cap, 1000 mL Capacity, 53 – 415 thread, 67 × 127 × 180 mm Dimension
Lot Number	VHC220321
Date of Manufacture (yyyy-mm-dd)	2021-03-22
Country of Origin	Manufactured in India
Date of Issue (yyyy-mm-dd)	2022-03-21

Quality System Compliance

Products are manufactured under the **ISO 9001:2015 & ISO 13485:2016** standard. Products are inspected and controlled through whole production processing in accordance with current applicable product specifications and QC SOP.

QC Testing

Representative products samples are collected and inspected in accordance with current applicable QC requirements.

Test Name	Test Reference	Accuracy	Result
Leak Test	QWI-QC-03	General Insp. Level-1, AQL -1.5	Leak Proof
Autoclaving	QWI-QC-07	General Insp. Level-1, AQL -1.5	NA
Dimensional	QSF-QCIR-02	General Insp. Level-1, AQL -1.5	Within Tolerance

Product Specifications

The product is in conformance with the requirements of the following standards and regulations:

Material **Carboy** High Density Polyethylene (HDPE).
Closure Polypropylene (PP).

Above material is 100% Virgin, Compliance with FDA 21 CFR

Temperature Resistance	up to +120 °C.
Sterilization	Not applicable.
ATP Assay	Not applicable.
DNase & RNase Free	Not applicable.
BSE/TSE	No use of any raw material produced from or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE).
Cytotoxicity	Not applicable.
Non-Pyrogenic Statement	Not applicable.
Latex Statement	Not applicable.
BPA Statement	Not applicable.
DEHP Statement	Please see RoHS information below.
RoHS	No substances (Lead, Cadmium, Mercury, Hexavalent Chromium (Cr6+), Poly Brominated Biphenyls (PBB), Poly Brominated Diphenyl ethers (PBDE), Bis(2-Ethyhexyl) Phthalate (DEHP), Benzyl Butyl Phthalate (BBP), Dibutyl Phthalate (DBP), Diisobutyl Phthalate (DIBP)) are used in manufacturing the raw materials and final product. No routinely analysis is performed.
REACH Statement	No substances as listed on >Substance of Very High Concern< candidate list published 2022-01-17 are used in manufacturing of the raw materials used to produce this product. No routinely analysis is performed.
Storage Conditions	Store in Room temperature.

Disclaimer: VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. The purpose of the product is for use in laboratory only.
 This document has been produced electronically and is valid without a signature