

Certificate of Conformity

VWR Catalogue Number	See List
Description	White Overalls with hood or collar
Country of Origin	Manufactured in China
Date of Issue (yyyy-mm-dd)	2021-03-12

Overalls with hood

VWR Catalogue Number	Size
113-1983	M
113-1984	L
113-1985	XL
113-1986	XXL
113-1987	3XL

Overalls with Collar

VWR Catalogue Number	Size
113-1979	L
113-1980	XL
113-1981	XXL
113-1982	3XL

Quality System Compliance

Products are manufactured under the [ISO 9001:2015 & ISO 13485:2003](#) standard. Products are inspected and controlled through whole production processing in accordance with current applicable product specifications and QC SOP.

QC Testing

Representative products samples are collected and inspected in accordance with current applicable QC requirements. Inspected according to AQL level II 4.0 sampling plan per MIL-STD-105E.

Inspected according to ASTM D5034 standard.

	MD	CD
Tensile Strength [N/5 cm]	> 85	> 60
Elongation rate [%]	>50	> 50

Product Specifications

Material	Spunbonded Polypropylene (SPP) 40 g/m ² .
Sterilization	Not applicable.
ATP Assay	Not applicable.
DNase & RNase Free	Not applicable.
BSE/TSE	No use of any raw material produced from or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE).
Cytotoxicity	Not applicable.
Non-Pyrogenic Statement	Not applicable.
Latex Statement	The above product is latex free.
BPA Statement	Not applicable.
DEHP Statement	Not applicable.
RoHS	Not applicable.
REACH Statement	Not applicable.
Storage Conditions	Room temperature, dry and shady.
Shelf Life	Not applicable.

Disclaimer: VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. The purpose of the product is for use in laboratory only.
This document has been produced electronically and is valid without a signature.