



Freudenberg Medical, LLC
1110 Mark Avenue
Carpinteria, CA 93013
Telephone: (805) 684-3304 Fax: (805) 684-1934

Date: 09/10/2025

Address:

Freudenberg Medical Europe
2-8 Liebigstrasse
67661 KAISERSLAUTERN
GERMANY

Freudenberg Medical
Material Number Batch Number

45639007 - 0000351297

Ware entspricht der Spezifikation
QS: 22. SEP. 2025

Subject: Certificate of Compliance

PO Number	: 4500041662	Quantity Shipped	: 6 EA
Delivery Number	: 80419177	Manufacturing Date	: 08/2025
Part Description	: HM Clamp Gasket, 0.75" ID (25 Pack)	Expiration date	: 08/2028

CUSTOMER P/N	Internal P/N	Internal REV	Internal L/N
45639007	70-400-07	A	0000351297
Internal D/N	Internal D/N Rev		
70-400	A		

Components

Material	Part Number	Lot Number
Manufacturer	Mfr Part Number	Mfr Lot Number
Dupont, Elastomer Silicone Q7-4780 A	30-011-01	0000324957 H033OBP000
Dupont, Elastomer Silicone Q7-4780 A	30-011-01	0000346251 H033P52018
Dupont, Elastomer Silicone Q7-4780 B	30-012-01	0000324958 H033OBP000
Dupont, Elastomer Silicone Q7-4780 B	30-012-01	0000346252 H033P52018
LBL Stock 3.75 X 7.5 Avery 2Sht	36-313-01	0000345979
HM Clamp Gasket, 0.75" ID	50-148-03	0000338930
HM Clamp Gasket, 0.75" ID	50-148-03	0000350525

Freudenberg Medical certifies that this product was manufactured in accordance with the material specifications and inspection criteria provided. The materials used in manufacture meet USP Class VI biocompatibility requirements. These platinum-cured silicone products have full traceability to raw materials and contain no peroxide catalysts.

This product is manufactured following Good Manufacturing Practices(GMP) in a FDA registered facility. Additionally, Freudenberg Medical has been certified to meet the requirements of ISO 13485.