

Disposable SW Plate Riplate® SW 96, PP, 2ml



Product description

Plate 96-SW (square well), 2 ml


Part Number	43001-0020
Well volume	2 ml
Well bottom	U-Bottom
Plate color	clear
Sterile	no
Re-usable	no
Patient contact	no

Materials

Well	PP – Polypropylene
-------------	--------------------

Dimensions

Plate Length	127.76 mm
Plate Height	44.00 mm
Plate Width	85.48 mm



Packaging

Box	20 plates
Pallet	112 boxes (2,240 plates)

Specifications are subject to change without prior notice.

In case of any questions, please contact our sales representative via phone, e-mail or our website.

<https://www.ritter-medical.de/en/contact>