


Physics

From classic and cutting-edge apparatus and demos, to time-saving activities, durable equipment, and essential supplies for your lab - we have what you need to turn physics lessons into physics connections.



Specifications:

Digital Multimeter E

DC voltage:	2 V – 600 V, 4 ranges, $\pm 0.8\%$ ± 5 digits
AC voltage:	2 V – 600 V, 4 ranges, $\pm 1\%$ ± 5 digits
DC current:	20 μA – 10 A, 7 ranges, $\pm 1.8\%$ ± 2 digits
AC current:	20 μA – 10 A, 7 ranges, $\pm 2\%$ ± 3 digits
Resistance:	200 Ω – 200 M Ω , 7 ranges, $\pm 1\%$ ± 4 digits
Display:	3½- display LCD, 24 mm, max. 1999
Operating voltage:	9 V battery 6F22
Fuse:	F1: F 2 A / 600 V F2: F 10 A / 600 V, $I_{\text{max.}} = 10$ A for 10 s with minimum interval 15 mins
Dimensions:	90x190x35 mm ³ approx.
Weight:	310 g approx. (with battery)