



Certificate of Conformity

VWR 10 EVENT DIGITAL THERMOMETER

VWR International certifies that the products comply in every respect with the relevant particular specifications and/or drawings and/or relevant standards and regulations in force, mentioned in this document.

VWR International certifie que les produits sont conformes à tous égards aux spécifications particulières appropriées et / ou dessins, normes et règlements en vigueur, mentionnées dans le présent document

VWR International bestätigt, dass das Produkte allen zutreffenden Spezifikations- und/oder Zeichnungs- und/oder Standardanforderungen sowie den im Dokument aufgeführten gültigen Vorschriften entspricht.

Catalog No.	Range (°C)	Probe	WxDxH (mm)	LxØ (mm)
620-2343	-50...+200 °C	Stainless steel	65x18.5x110	125x3.5

Accuracy (°C): ±0,5 (-20...+100); ±1,5 (<+100) °C

Resolution (°C): 0.1

Other Characteristics:

- Date and time of up to ten events that fall outside of the limits.
- Audible alarm, low or high alarm icon, indicates maximum/minimum temperatures.
- Verification model with glass bead filled bottle to eliminate detection of rapid changes and protect from contamination.
- Include 10 event recording, minimum/maximum temperature readings, timers and alarms.

Disclaimer: VWR wants to state that this declaration will not discharge the user from his obligation to ensure the product is suitable for the intended use. - If there is any doubt ask us.

Date of issue: **June 2017**