

A



B



C



D



MICRO TIPS

Many of our low volume tips are available with our MicroPoint™ design. This design minimizes surface area around the tip's orifice, reducing surface tension and allowing the most accurate dispense possible - a must when dispensing less than 2 μ l.

VOLUME 10 μ l STANDARD ART TIPS

Catalog Number	Product	Description	Unit Packaging
2139	ART 10	10 μ l w/ART barrier MicroPoint Design	Racked, pre-sterilized, 96 tips/tray, 10 trays/pack
2139B	ART 10	10 μ l w/ART barrier MicroPoint Design	Bulk, non-sterile, 1000 tips/bag
2139-RI	ART 10	10 μ l Pipette Tip Reload Insert w/ART barrier MicroPoint Design	Reload, pre-sterilized, 96 tips/tray, 10 trays/pack
2139-HR	ART 10	10 μ l Pipette Tip Hinged Rack w/ART barrier MicroPoint Design	Hinged Racked, pre-sterilized, 96 tips/tray, 10 trays/pack