

Certificate of Conformity

VWR Catalogue Number	See List
Description	Visitor coats, Velcro® closure, VWR® Basic with Weight 30 g/m ²
Country of Origin	Manufactured in China
Date of Issue (yyyy-mm-dd)	2022-07-04

Catalogue Number	Colour	Clothing Size	Length [cm]
113-8955	White	XS	92
113-8956	White	S	96
113-8957	White	M	100
113-8958	White	L	104
113-8959	White	XL	108
113-8960	White	XXL	112
113-8962	White	3XL	116
113-8963	Blue	XS	92
113-8964	Blue	S	96
113-8965	Blue	M	100
113-8966	Blue	L	104
113-8967	Blue	XL	108
113-8969	Blue	XXL	112
113-8971	Blue	3XL	116

Quality System Compliance

Products are manufactured under the **ISO 9001:2015 & ISO 13485:2016** standard. Products are inspected and controlled through whole production processing in accordance with current applicable product specifications and QC SOP.

QC Testing

Representative products samples are collected and inspected in accordance with current applicable QC requirements. Inspected according to AQL level II 4.0 sampling plan per MIL-STD-105E.

	Test Method	MD	CD
Break Tensile Properties [N/5cm]	ISO13934-1	> 50	> 30
Elongation rate [%]	ISO13934-1	> 30	> 30

Product Specifications

Material	30gsm SPP-Blue (Spunbonded polypropylene fabric). The raw materials used to produce the products are listed in List of Authorised Substances of Regulation EU/10/2011. The material grades conform to EU Directive 1935/2004 on materials and articles intended to come into contact with food.
Sterilization	Not Applicable.
ATP Assay	Not Applicable.
DNase & RNase Free	Not Applicable.
BSE/TSE	No use of any raw material produced from or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE).
Cytotoxicity	Not Applicable.
Non-Pyrogenic Statement	Not Applicable.
Latex Statement	This product is not made from natural rubber latex.
BPA Statement	Not Applicable.
DEHP Statement	Not Applicable.

RoHS	Not Applicable.
REACH Statement	Not Applicable.
Storage Conditions	Room temperature, dry and shady.
Shelf Life	Not Applicable.

Disclaimer: VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. The purpose of the product is for use in laboratory only.
This document has been produced electronically and is valid without a signature.