

Material	23835.294
Material description	Petroleum spirit 40-60°C (max. 0.01% aromatic hydrocarbons)
Grade	AnalaR NORMAPUR ACS
CAS Number	64742-49-0

Characteristics	Specifications
Acidity	≤ 0.0002 meq/g
Bromine value	≤ 1.0
Colouration	≤ 10 APHA
Boiling range (40-60°C) (Vol %)	≥ 90 %
Substances coloured by H ₂ SO ₄	≤ 150 APHA
Benzene	≤ 100 ppm
Evaporation residue	≤ 10 ppm
Total S (as SO ₄)	≤ 50 ppm
Water	≤ 100 ppm
Al (Aluminium)	≤ 0.2 ppm
Ba (Barium)	≤ 0.1 ppm
Ca (Calcium)	≤ 0.5 ppm
Cd (Cadmium)	≤ 0.05 ppm
Co (Cobalt)	≤ 0.02 ppm
Cr (Chromium)	≤ 0.02 ppm
Cu (Copper)	≤ 0.02 ppm
Fe (Iron)	≤ 0.1 ppm
K (Potassium)	≤ 0.3 ppm
Mg (Magnesium)	≤ 0.1 ppm
Mn (Manganese)	≤ 0.02 ppm
Na (Sodium)	≤ 1 ppm
Ni (Nickel)	≤ 0.02 ppm
Pb (Lead)	≤ 0.1 ppm
Sn (Tin)	≤ 0.1 ppm
Sr (Strontium)	≤ 0.03 ppm
Zn (Zinc)	≤ 0.1 ppm
Conforms to ACS	Passes test



Product Specification

Signature

We certify that this batch conforms to the specifications listed above.

This document has been produced electronically and is valid without a signature.

Isabelle Habay, Head of Laboratory - Briare
VWR International S.A.S.; Z.I. de Vaugereau; FR-45250
Briare; France