

Certificate of Conformity

VWR Catalogue Number	See List
Description	PCR plates, 96-well with Coloured Lettering
Country of Origin	Manufactured in Germany
Date of Issue (yyyy-mm-dd)	2021-09-27

Catalogue Number	Description	Well Colour	Well Volume [mL]
732-2387	Standard profile, non skirted	Clear	0.20
732-2390	Standard profile, half skirted	Clear	0.20
732-2386	Low profile, non skirted	Clear	0.15
732-2388	Low profile, half skirted	Clear	0.15
732-2389	Low profile, raised half skirted	Clear	0.15
732-3388	Standard profile, non skirted	White	0.20
732-3387	Low profile, non skirted	White	0.15
732-3389	Low profile, half skirted	White	0.15
732-3390	Low profile, raised half skirt	White	0.15

Quality System Compliance

Products are manufactured under the **ISO 9001:2015** standard. Products are inspected and controlled through whole production processing in accordance with current applicable product specifications QC SOP and Clean room Class 8 according to ISO 14644-1.

QC Testing

Representative products samples are collected and inspected in accordance with current applicable QC requirements.

Product Specifications

Material	Polypropylene.
Sterilization	Not Applicable.
ATP Assay	Not Applicable.


Endotoxin Free	This product is free of Endotoxin.
DNase & RNase Free	This product is free of any detectable DNase/RNase contamination.
PCR Inhibitors	This product is free of PCR Inhibitors.
BSE/TSE	Tallow derived additives may be used in the manufacture of this product. "Mad Cow" tallow derived materials used in this product fulfill the requirements laid down in the Regulations 1069/2009/EC, and 142/2011/EC, and the "Note for Guidance EMA/410/01, and as amended.
Cytotoxicity	Not Applicable.
Non-Pyrogenic Statement	The acceptance level for product is 0.01 EU/mL device (LAL Gel Clot Method).
Latex Statement	This product is not made from natural rubber latex. Bisphenol are not used in the manufacture of the raw material and are not expected to be present.
DEHP Statement	Please see RoHS information below.
RoHS	No substances (Lead, Cadmium, Mercury, Hexavalent Chromium (Cr6+), Poly Brominated Biphenyls (PBB), Poly Brominated Diphenyl ethers (PBDE), Bis(2-Ethyhexyl) Phthalate (DEHP), Benzyl Butyl Phthalate (BBP), Dibutyl Phthalate (DBP), Diisobutyl Phthalate (DIBP)) are used in manufacturing the raw materials and final product. No routinely analyse is performed.
REACH Statement	No substances as listed on >Substance of Very High Concern< candidate list published 2021-07-08 are used in manufacturing of the raw materials used to produce this product. No routinely analysis is performed.
Storage Conditions	Store in a Room temperature.
Shelf Life	5 years.

Disclaimer: VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. The purpose of the product is for use in laboratory only.
 This document has been produced electronically and is valid without a signature.

Annex


Technical Drawing

732-2387



L	120.0 mm	Yn	71.5 mm
W	80.0 mm	P	9.0 mm
H	20.7 mm	D	5.4 mm
X1	97.5 mm		
Xn	109.5 mm	C	12
Y1	63.5 mm	R	8


732-2390



L	127.8 mm
W	85.5 mm
H	20.7 mm
X1	14.4 mm
Xn	113.4 mm
Y1	11.2 mm


Yn	74.2 mm
P	9.0 mm
D	5.4 mm
C	12
R	8

732-2386




L	120.0 mm	Y _n	71.5 mm
W	80.0 mm	P	9.0 mm
H	15.3 mm	D	5.2 mm
X ₁	10.5 mm		
X _n	109.5 mm	C	12
Y ₁	8.5 mm	R	8

732-2388




L	127.8 mm	Yn	74.2 mm
W	85.5 mm	P	9.0 mm
H	15.3 mm	D	5.4 mm
X1	14.4 mm		
Xn	114.4 mm	C	12
Y1	11.2 mm	R	8

732-2389




L	127.2 mm	Yn	74.1 mm
W	85.2 mm	P	9.0 mm
H	15.3 mm	D	5.4 mm
X1	14.1 mm		
Xn	113.1 mm	C	12
Y1	11.1 mm	R	8

732-3388




L	120.0 mm	Y_n	71.5 mm
W	80.0 mm	P	9.0 mm
H	20.7 mm	D	5.4 mm
X_1	97.5 mm		
X_n	109.5 mm	C	12
Y_1	63.5 mm	R	8

732-3387




L	120.0 mm	Yn	71.5 mm
W	80.0 mm	P	9.0 mm
H	15.3 mm	D	5.4 mm
X1	10.5 mm		
Xn	109.5 mm	C	12
Y1	8.5 mm	R	8

732-3389



L	127.8 mm	Yn	74.2 mm
W	85.5 mm	P	9.0 mm
H	15.3 mm	D	5.4 mm
X1	14.4 mm		
Xn	114.4 mm	C	12
Y1	11.2 mm	R	8

732-3390



L	127.2 mm	Yn	74.1 mm
W	85.2 mm	P	9.0 mm
H	15.3 mm	D	5.4 mm
X1	14.1 mm		
Xn	113.1 mm	C	12
Y1	11.1 mm	R	8