

	TECHNICAL DATA SHEET	300200EN	Fecha: 10/12/13
		Revisión 7	Página 1/2

Codes: 300200, 310200, 300201, 300202, 310202, 300203, 310253.1
Sterile swabs, individually wrapped

Page 18 in 2014 General Catalogue

Code	Description	Package	Case quantity	Case weight	Case volume
300200	Wood + cotton	Peel-pack	2x1000	2.6	0.029
310200	Wood + cotton	Flow-pack	2x1000	2.6	0.029
300201	Snappable PS + cotton	Peel-pack	2x1000	3.5	0.029
300202	Snappable PS + viscose	Peel-pack	2x1000	3.5	0.029
310202	Snappable PS + viscose	Flow-pack	2x1000	2.02	0.029
300203	Aluminium + cotton	Peel-pack	2x1000	2.4	0.029
310253.1	Aluminium + viscose	Flow-pack	2x1000	2.3	0.029

Directive 93/42/EEC. Medical devices. Class IIa.

SPECIFIC USE AND DESCRIPTION

Sterile swabs for clinical use, designed for the collection of biological samples. Mainly used for the preparation of smears (stained and then observed through a microscope). Also useful for isolations in culture media.

Expiry date: 4 years from sterilization date.

Certificates of sterilization are available on request. The sterilization method is printed on the packaging.

Swabs are also useful for surface spreading techniques, required in antibiotic susceptibility testing for instance.

Other applications include the cleaning of biological cavities or places hard to reach.

Swabs for the collection of microbiological samples in open wounds are ISO 9001 manufactured and have the CE Marking 0318, Class IIa. Copies of both certificates are available upon request.

Deltalab is one of the few European firms to hold the CE Mark Class IIa, being therefore allowed to sell its swabs to other EC members and to those countries demanding the highest product quality.



Created	Verified	Approved
		

	TECHNICAL DATA SHEET	300200EN	Fecha: 10/12/13
		Revisión 7	Página 2/2

PACKAGING DESCRIPTION

- **Peel-pack:** bag of 1 unit. Fibre-free peeling paper. One side made of paper; the other one, made of two films: one of PVC and another one of polyethylene.
- **Flow-pack:** bioriented polypropylene bag of 1 unit. A low cost and safe packaging.

On both packagings types the following data are printed: product code, description and brand, batch number and expiry date.

Code	Shaft + tip	Shaft length mm	Shaft diameter mm	Tip diameter mm	Length shaft+tip mm
300200 (peel) 310200 (flow)	Wood + cotton	150	2.5	5.0	150 (136+14)
300201 (peel)	Snappable PS + cotton	150	2.5	5.0	150 (136+14)
300202 (peel) 310202 (flow)	Snappable PS + viscose	150	2.5	5.0	150 (136+14)
300203 (peel)	Aluminium + cotton	150	0.9	1.5	150 (137+13)
310253.1 (flow)	Aluminium + viscose	150	0.9	1.5	150 (137+13)

Created	Verified	Approved
		