

# Certificate of Conformity

|                                               |                                                                                              |
|-----------------------------------------------|----------------------------------------------------------------------------------------------|
| <b>VWR Catalogue Number</b>                   | 212-0726                                                                                     |
| <b>Description</b>                            | Test tubes with red cap, sterile,<br>5 ml Capacity, 75 mm Length, 13 mm<br>External Diameter |
| <b>Lot Number</b>                             | 21F1773                                                                                      |
| <b>Sterilization Expiry Date (yyyy-mm-dd)</b> | 2027-08                                                                                      |
| <b>Country of Origin</b>                      | Manufactured in Italy                                                                        |
| <b>Date of Issue (yyyy-mm-dd)</b>             | 2022-11-07                                                                                   |

## Quality System Compliance

Products are manufactured under the **ISO 9001:2015** standard. Products are inspected and controlled through whole production processing in accordance with current applicable product specifications and QC SOP.

## QC Testing

Representative products samples are collected and inspected in accordance with current applicable QC requirements.

## Product Specifications

The product is in conformance with the requirements of the following standards and regulations:

|                               |                                                                                                                        |
|-------------------------------|------------------------------------------------------------------------------------------------------------------------|
| <b>Material</b>               | Polystyrene.<br><br>The products are considered safe to use in food, pharmaceutical and cosmetic contact applications. |
| <b>Sterilization</b>          | Products are sterilized using Beta rays and meets SAL 10 <sup>-6</sup> with a specified dosage range of 17 – 20 kG.    |
| <b>ATP Assay</b>              | Not applicable.                                                                                                        |
| <b>DNase &amp; RNase Free</b> | Not applicable.                                                                                                        |

|                                |                                                                                                                                                                                                                                                                                                                                                             |
|--------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>BSE/TSE</b>                 | No use of any raw material produced from or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE). |
| <b>GMO Statement</b>           | Genetically modified materials (GMO) have not been intentionally incorporated in the manufacturing of the product.                                                                                                                                                                                                                                          |
| <b>Cytotoxicity</b>            | Not applicable.                                                                                                                                                                                                                                                                                                                                             |
| <b>Non-Pyrogenic Statement</b> | Not applicable.                                                                                                                                                                                                                                                                                                                                             |
| <b>Latex Statement</b>         | The products are not made from natural rubber latex.                                                                                                                                                                                                                                                                                                        |
| <b>BPA Statement</b>           | Bisphenol A (BPA) has not been intentionally incorporated in the manufacturing of the products.                                                                                                                                                                                                                                                             |
| <b>DEHP Statement</b>          | DEHP has not been intentionally incorporated in the manufacturing of the products.                                                                                                                                                                                                                                                                          |
| <b>RoHS</b>                    | Not applicable.                                                                                                                                                                                                                                                                                                                                             |
| <b>REACH Statement</b>         | No substances as listed on >Substance of Very High Concern< candidate list published 2022-06-10 are used in manufacturing of the raw materials used to produce this product. No routinely analysis is performed.                                                                                                                                            |
| <b>Storage Conditions</b>      | Store in a dry place, away from sources of heat and direct sunlight.                                                                                                                                                                                                                                                                                        |

Disclaimer: VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. The purpose of the product is for use in laboratory only.  
This document has been produced electronically and is valid without a signature