



Certificate of Conformity

VWR DIGITAL THERMOMETERS

VWR International certifies that the products comply in every respect with the relevant particular specifications and/or drawings and/or relevant standards and regulations in force, mentioned in this document.

VWR International certifie que les produits sont conformes à tous égards aux spécifications particulières appropriées et / ou dessins, normes et règlements en vigueur, mentionnées dans le présent document

VWR International bestätigt, dass das Produkte allen zutreffenden Spezifikations- und/oder Zeichnungs- und/oder Standardanforderungen sowie den im Dokument aufgeführten gültigen Vorschriften entspricht.

Catalog No.	Range (°C)	Probe	Pk
620-0904	-50...+200 °C	Glass bead recyclable bottle	1
620-0916	-50...+200 °C	Removable stainless steel	1

Accuracy (°C): ±0,5 (-20...+70); ±1,0 (70...120); ±2,0 (-50...-20; 120...150); ±4,0 (150...200) °C °C

Resolution (°C): 0.1

Weight (g) : 136

Other Characteristics:

- These thermometers allow for measurements in food and beverage processing and storage.
- LCD display.
- 'HOLD' function.
- Maximum/minimum temperature.

Disclaimer: VWR wants to state that this declaration will not discharge the user from his obligation to ensure the product is suitable for the intended use. - If there is any doubt ask us.

Date of issue: **March 2017**