



Specification

1.03979.0500 Meat extract dry, granulated, for microbiology

Specification		
Appearance	Lightly brownishyellow granulate	
pH-value (5 %; water)	6.0 - 7.5	
Total nitrogen (N) (Kjeldahl)	11.5 - 12.5	%
Sulfated ash (800 °C)	≤ 18.0	%
Loss on Drying (105 °C)	≤ 6.0	%
Amino nitrogen (as N)	3.5 - 4.5	%
Nitrite (NO ₂)	passes test	
Suitability for microbiology	passes test	
Identity (NIR)	passes test	

Stefanie Fischer

Responsible laboratory manager quality control

This document has been produced electronically and is valid without a signature.