

VWR Brand Autoclavable Biohazard Bags

VWR CATALOG NUMBERS:

[Manufacturing country of origin for all below: USA, HS code 392329]

2 mil (50.8μ) PP Autoclavable Biohazard Bags				size	
Red Bags	Clear Bags	Orange Bags	Yellow Bags	mm	inches
VWRI 129-1550		VWRI 129-1564	VWRI 129-1148	200 x 305	8 x 12
VWRI 129-1551		VWRI 129-0024	VWRI 129-1149	355 x 480	14 x 19
	VWRI 129-0028			305 x 610	12 x 24
VWRI 129-1552		VWRI 129-0025	VWRI 129-1150	480 x 610	19 x 23
	VWRI 129-1574				24 x 30
	VWRI 129-0029		VWRI 129-1151		24 x 36
VWRI 129-1553		VWRI 129-0026		610 x 760	25 x 35
VWRI 129-1554		VWRI 129-0027		635 x 889	35 x 45
	VWRI 129-0029		VWRI 129-1148	950 x 1220	37 x 48

DESCRIPTION:

VWR Autoclavable Biohazard bags made from Suprop® Polypropylene film, specifically engineered to exhibit seal ability, exceptional toughness and excellent autoclavable temperature resistance (typical program of 20 mins. @ Operating temperature to 135°C, film will melt above 149°C) as utilized in some packaging and sterilization applications.

TYPICAL PRODUCT CHARACTERISTICS:

Values expressed below for available nominal film thicknesses.

PROPERTIES		2.0 mil (50.8μ)	Units	Test Method
Elmendorf Tear Resistance *	MD ↕	200	grams	ASTM D1922
	TD ↔	980	grams	
Dart Impact Strength		>400	grams	ASTM D1709
Coefficient of Friction		<0.20		ASTM D1894
Yield		15300	sq.in./lb.	Spartech-Tufpak

FOR YOUR PROTECTION:

The information and recommendations in this publication are, to the best of our knowledge, accurate and reliable. Suggestions made concerning uses or applications are only the opinion of Spartech-Tufpak, and users should make their own tests to determine the suitability of these products for their own purposes. However, because of numerous factors affecting results, Spartech-Tufpak. makes no warranty of any kind, express or implied, including those of merchantability and fitness for purpose, other than that the material conforms to its Standard Specifications. Statements herein, therefore, should not be construed as representations or warranties. The responsibility of Spartech-Tufpak. for claims arising out of breach of warranty, negligence, strict liability, or otherwise is limited to the purchase price of the material. Statements concerning the use of products described herein are not to be construed as recommending the infringement of any patent and no liability for infringement arising out of any such use is assumed. 12/2021