

Contec® Products for Mopping Cleanroom Environments

Contec® Cleanroom has a wide range of mopping products for all areas of the cleanroom. This document contains the ideal products for critical environments because when used together:

- Sealed edge flat foam mop provides even disinfectant application to surfaces
- Layflat hardware makes it easy to clean under and around equipment
- Lightweight handles allow for easy cleaning of ceilings and other hard-to-reach areas
- Hands-free mop removal to prevent cross contamination

VertiKlean® MAX™ Sealed Edge Mop

- Laser sealed reducing the chance of fiber and particle contamination
- Less than 2 inch clearance under and behind equipment
- Compatible with a wide range of disinfectants and solutions
- Hydrophilic polyurethane foam retains liquids to evenly release without dripping



Part No.	Description		Size	Packaging
10803-094	Medium VertiKlean MAX Sterile sealed edge mop head	 Sterile	15" x 5" x 0.375" (38 x 13 x 1 cm)	1/bag; 24 bags/case
76199-908	Large VertiKlean MAX Sterile sealed edge mop head	 Sterile	18" x 5" x 0.375" (46 x 13 x 1 cm)	1/bag; 20 bags/case
10171-158	Medium VertiKlean MAX Sealed edge mop head		15" x 5" x 0.375" (38 x 13 x 1 cm)	4/bag; 4 bags/case
76199-898	Large VertiKlean MAX Sealed edge mop head		18" x 5" x 0.375" (46 x 13 x 1 cm)	4/bag; 5 bags/case
76299-418	QuickConnect™ Layflat Mop Frame		11.5" (29 cm)	1/case



Order on
avantorsciences.com

QuickConnect™ Carbon Fiber Handles

- Fully autoclavable
- Compatible with all Contec mops
- Compatible with common cleanroom disinfectants and solutions
- Multiple lengths from 18 inches to 28'
- Lightweight compared to stainless steel



Part No.	Description	Size	Packaging
77322-296	QuickConnect Carbon Fiber Telescopic Handle	40" - 70" (102 - 178 cm)	1/case
77322-288	QuickConnect Carbon Fiber Telescopic Handle	63" - 114" (160 - 290 cm)	1/case
77322-290	QuickConnect Carbon Fiber Telescopic Handle	66" - 224" (168 - 568 cm)	1/case
77322-294	QuickConnect Carbon Fiber Telescopic Handle	69" - 334" (175 - 848 cm)	1/case

Note: QuickConnect handles are also available in stainless steel, anodized aluminum and fiberglass.

Compact Bucket System

- Fully autoclavable
- Hands-free removal tool to prevent operators from touching used mop heads therefore preventing cross contamination
- Storage for extra supplies including mops, wipes and handles
- Quiet, lockable casters move easily



Part No.	Description	Size	Packaging
77786-112	Two 15 L polypropylene buckets, with VertiKlean MAX mop removal tool attached, utility handle, two 6L buckets and one debris screen	25" x 26.5" x 41" (63 x 67 x 104 cm)	1/case
77786-114	Three 15L polypropylene buckets, with VertiKlean MAX mop removal tool attached, utility handle, two 6L buckets and one debris screen	32" x 26.5" x 41" (81 x 67 x 104 cm)	1/case
76199-932	Stainless Steel Flat Mop Wringer	-	1/case

avantor
ScienceCentral™
avantorsciences.com

Any prices, product, and/or services details are current when published and subject to change. Offers subject to product availability. Void where prohibited by law or company policy. Any claims concerning products manufactured by third parties are the sole claims of the manufacturer and not those of Avantor or its affiliates. All products are subject to the terms and conditions of sale of Avantor or its affiliates, including the applicable limitations of liability contained therein. The information herein should not be considered a warranty of any kind. Please carefully review any product information sheets and safety information for the product. Trademarks are owned by Avantor Inc. or its affiliates, unless otherwise noted. | Visit [avantorsciences.com](https://www.avantorsciences.com) to view our privacy policy and additional disclaimers. © 2025 Avantor Inc. All rights reserved.