

Certificate of Conformity

VWR Catalogue Number	See List
Description	Disposable Inoculating loop
Country of Origin	Manufactured in Italy
Date of Issue (yyyy-mm-dd)	2023-09-07

Catalogue Number	Type	Capacity [µl]	Colour
612-9351	Loops, rigid	1	Dark Green
612-9352	Loops, rigid	1	Dark Green
612-9353	Loops, rigid	10	Dark Blue
612-9354	Loops, rigid	10	Dark Blue
612-9355	Loops, flexible	1	Light Green
612-9356	Loops, flexible	1	Light Green
612-9357	Loops, flexible	10	Light Blue
612-9358	Loops, flexible	10	Light Blue

Quality System Compliance

Products are manufactured under the **ISO 9001:2015 & ISO 13485:2016** standard. Products are inspected and controlled through whole production processing in accordance with current applicable product specifications and QC SOP.

QC Testing

Representative products samples are collected and inspected in accordance with current applicable QC requirements.

Product Specifications

Material Polystyrene.

Sterilization Product labelled as sterile is gamma irradiated. Sterilisation procedure is validated in conformity with BS EN ISO 11137-1:2015 and BS EN ISO 11137-2:2015 to guarantee SAL 10⁻⁶.

ATP Assay	Not applicable.
DNase & RNase Free	Not applicable.
BSE/TSE	No use of any raw material produced from or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE).
Cytotoxicity	Not applicable.
Non-Pyrogenic Statement	Not applicable.
Latex Statement	This product is not made from natural rubber latex.
BPA Statement	Not applicable.
DEHP Statement	Not applicable.
RoHS	Not applicable.
REACH Statement	No substances as listed on >Substance of Very High Concern< candidate list published 2023-06-14 are used in manufacturing of the raw materials used to produce this product. No routinely analysis is performed.
Storage Conditions	No special storage condition required.
Shelf Life	3 years.

Disclaimer: VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. The purpose of the product is for use in laboratory only.
This document has been produced electronically and is valid without a signature