



## Certificate of Analysis

**Product:** 0.9% SODIUM CHLORIDE IRRIGATING SOLUTION

**Lot #: G123166**  
**Code: 2F7124**  
**Expiry: MAY 2019**

**Manufactured Date: 05/09/2016**  
**Sterilization Date: 05/09/2016**  
**Mixing Specification: 03-15-19-156**  
**Issue Date: 07/02/2002**

| TEST NUMBER                                                 | TEST NAME                                       | LIMIT                            | RESULT         |
|-------------------------------------------------------------|-------------------------------------------------|----------------------------------|----------------|
| 11-07-19-001                                                | NaCl (g/L)                                      | 8.55 - 9.45 g/L                  | 9.05 g/L       |
| 11-20-19-006                                                | Sodium ID- USP                                  | Positive                         | Positive       |
| 11-20-19-007                                                | Sodium ID- Flame                                | Positive                         | Positive       |
| 11-21-16-003                                                | pH at 25 deg. C                                 | 4.5 - 7.0                        | 5.7            |
| 20-03-07-006                                                | Particle Analysis - Large Volume (Instrumental) |                                  |                |
|                                                             | Particles >= 10 micron                          | NMT 25 particles/mL >=10 microns | 0 ct/mL        |
|                                                             | Particles >= 25 micron                          | NMT 3 particles/mL >=25 microns  | 0 ct/mL        |
| 13-01-10-313                                                | Bacterial Endotoxin                             | 0.5 EU/mL                        | < 0.0100 EU/mL |
| 14-04-01-120                                                | Sterility                                       | Pass Parametric Release          | Pass           |
| PARAMETRIC RELEASE IN ACCORDANCE WITH U.S. FDA REQUIREMENTS |                                                 |                                  |                |

THIS PRODUCT WAS MANUFACTURED ACCORDING TO CURRENT GMP,  
FULLFILLS THE REQUIREMENTS OF THE MASTER PRODUCTION DOCUMENT AND TESTED TO USP COMPLIANCE

E. Soliz 5-11-16  
Prepared By/Date

Ellen Soliz

Twila Proctor 5/11/16  
Quality Management Approval /Date  
Twila Proctor