

1. Part Number: 1095S-150

2. Photo:



3. Product Information:

Description	V-Spatula 145mm x 25mm
Method of Construction:	Crevice free, all welds ground and polished
Materials:	316L stainless steel
Surface Finish:	Better than 0.5 microns Ra
Overall Length:	227mm
Blade Length:	145mm
Maximum Blade Width:	25mm
Nominal Weight:	160g
Recommended Storage Conditions:	Dry and ambient temperature
First Use:	This sampler should be thoroughly cleaned before being used for the first time

4. Approval:

Signed:



Name: Kuldeep Flora
Position: Quality Manager
Date: 26-Nov-19