

# Certificate of Conformity

|                                         |                                                                                                              |
|-----------------------------------------|--------------------------------------------------------------------------------------------------------------|
| <b>VWR Catalogue Number</b>             | 113-8226                                                                                                     |
| <b>Description</b>                      | Coveralls with Hood, with Elasticated cuffs, Ankles and waist, Thumb Loop, White Colour, Size M, Non-Sterile |
| <b>Lot Number</b>                       | 200910                                                                                                       |
| <b>Date of Manufacture (yyyy-mm-dd)</b> | 2020-09-10                                                                                                   |
| <b>Country of Origin</b>                | Manufactured in China                                                                                        |
| <b>Date of Issue (yyyy-mm-dd)</b>       | 2026-03-25                                                                                                   |

## Quality System Compliance

This certificate confirms that the product has been sourced from a supplier or manufacturer who has declared their Quality Management System (QMS) complies with **ISO 9001:2015 & ISO 13485:2016**.

## QC Testing

Representative products samples are collected and inspected in accordance with current applicable QC requirements. Inspected according to AQL level II 4.0 sampling plan per MIL-STD-105E.

Inspected according to ASTM D5034 standard.

|                                | MD  | CD |
|--------------------------------|-----|----|
| Break Tensile Strength [N]     | 116 | 77 |
| Elongation rate [%]            | 36  | 48 |
| Trapezoid Tearing Strength (N) | 118 | 86 |

## Product Specifications

|                               |                                                                                       |
|-------------------------------|---------------------------------------------------------------------------------------|
| <b>Material</b>               | 60g Spunbonded Polypropylene fabric/Meltblown /Spunbonded Polypropylene fabric (SMS). |
| <b>Sterilization</b>          | Non-Sterile.                                                                          |
| <b>ATP Assay</b>              | Not Applicable.                                                                       |
| <b>DNase &amp; RNase Free</b> | Not Applicable.                                                                       |

|                                |                                                                                                                                                                                                                                                                                                                                                             |
|--------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>BSE/TSE</b>                 | No use of any raw material produced from or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE). |
| <b>Cytotoxicity</b>            | Not Applicable.                                                                                                                                                                                                                                                                                                                                             |
| <b>Non-Pyrogenic Statement</b> | Not Applicable.                                                                                                                                                                                                                                                                                                                                             |
| <b>Latex Statement</b>         | This product is not made from natural rubber latex.                                                                                                                                                                                                                                                                                                         |
| <b>BPA Statement</b>           | Not Applicable.                                                                                                                                                                                                                                                                                                                                             |
| <b>DEHP Statement</b>          | Not Applicable.                                                                                                                                                                                                                                                                                                                                             |
| <b>RoHS</b>                    | Not Applicable.                                                                                                                                                                                                                                                                                                                                             |
| <b>REACH Statement</b>         | Not Applicable.                                                                                                                                                                                                                                                                                                                                             |
| <b>Storage Conditions</b>      | Store at 0° to 50°C, 45% to 65%.                                                                                                                                                                                                                                                                                                                            |
| <b>Shelf Life</b>              | Not Applicable.                                                                                                                                                                                                                                                                                                                                             |

Disclaimer: VWR's Quality Management System (QMS) is ISO 9001:2015 certified as a wholesale distributor of finished products, ensuring effective oversight of finished product handling, traceability, and distribution processes. VWR is not a manufacturer and does not engage in manufacturing, co-manufacturing or processing activities such as design, packaging, labeling, testing or quality release. Consequently, the supplier or manufacturer's QMS certification may differ from VWR's.

VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. This product is intended for research use only. It is not for use in diagnostic procedures, therapeutic applications, or any human or clinical use.

This document has been produced electronically and is valid without a signature.