



Product Specification Sheet

VWR SPECIMEN CONTAINER 120 ml, PP, FOR VACUUM TUBE, STERILE

VWR International certifies that the specified products meet all VWR specifications and quality criteria for release. The products comply in every respect with the relevant particular specifications, drawings, relevant standards and regulations in force.

VWR International certifie que les produits mentionnés respectent toutes les spécifications et critères de qualité VWR. Les produits satisfont aux spécifications particulières, aux plans, aux normes et règlements en vigueur s'y rapportant.

VWR International bestätigt hiermit, daß die angegebenen Produkte alle VWR Spezifikationen und Qualitätskriterien zur Freigabe erfüllen. Die Produkte entsprechen in jeder Hinsicht den maßgeblichen Spezifikationen, den Zeichnungen sowie den relevanten Normen und gültigen Vorschriften.

VWR Catalog No: **216-1161**

Cup Material: **Polypropylene**

Cap Material: **Polyethylene**

Appearance: **Cylindrical**

Outside sizes: $\varnothing$ lower: **47.9 mm**
 $\varnothing$ higher: **56 mm**

Height without cap: **71.5 mm**

Height with cap: **71.8 mm**

Colour: Cup: **Neutral, transparent**
Cap: **Yellow**

Other Characteristics: **Graduated at: 10, 20, 30, 40, 50, 60, 70, 80, 90, 100 ml**

Sterilization: **Sterile, Individually Packed**

Box Quantity: **250 pcs box**

Box Volume: **39 x 59 x 39 cm**

Box Weight: **4.60 kg**

Storage: **Do not expose or leave near open flames, do not leave the product around as it is not biodegradable.**

Disclaimer: VWR wants to state that these specifications can only be regarded as general information and will not discharge the user from his obligation to ensure the product is suitable for the intended use - If there is any doubt ask us.

Date of issue: **April 2013**