



Certificate of Conformity

VWR LABORATORY JACKS

VWR International certifies that the products comply in every respect with the relevant particular specifications and/or drawings and/or relevant standards and regulations in force, mentioned in this document.

VWR International certifie que les produits sont conformes à tous égards aux spécifications particulières appropriées et / ou dessins, normes et règlements en vigueur, mentionnées dans le présent document

VWR International bestätigt, dass das Produkte allen zutreffenden Spezifikations- und/oder Zeichnungs- und/oder Standardanforderungen sowie den im Dokument aufgeführten gültigen Vorschriften entspricht.

Catalog No.	Plate size (mm)	Max. static load (kg) *	Working height (mm)	Colour	Pk
245-0004	50x40	20	31 - 93	Green	1
245-0005	80x75	30	49 - 147	Blue	1
245-0006	120x140	50	60 - 243	Pink	1
245-0007	122x150	80	80 - 278	Yellow	1
245-0008	180x216	90	95 - 398	Blue	1
245-0009	240x320	90	120 - 500	Grey	1

* **Static load:** Maximum load once the height is adjusted

Material : Anodised aluminium

Other Characteristics:

- Smooth adjustment of working height.
- Extreme stability.
- Exceptionally sturdy.
- Corrosion resistant and non susceptible to a wide range of chemicals.

Disclaimer: VWR wants to state that this declaration will not discharge the user from his obligation to ensure the product is suitable for the intended use. - If there is any doubt ask us.

Date of issue: **June 2016**