



Product Specification Sheet

VWR GLASS FIBRE FILTERS, GRADE 693

VWR International certifies that the specified products meet all VWR specifications and quality criteria for release. The products comply in every respect with the relevant particular specifications, drawings, relevant standards and regulations in force.

VWR International certifie que les produits mentionnés respectent toutes les spécifications et critères de qualité VWR. Les produits satisfont aux spécifications particulières, aux plans, aux normes et règlements en vigueur s'y rapportant.

VWR International bestätigt hiermit, daß die angegebenen Produkte alle VWR Spezifikationen und Qualitätskriterien zur Freigabe erfüllen. Die Produkte entsprechen in jeder Hinsicht den maßgeblichen Spezifikationen, den Zeichnungen sowie den relevanten Normen und gültigen Vorschriften.

Catalogue N°: **516-0872**

Intended uses:

- Suspended solids analysis in wasted water
- Filtration of proteins, drinking water, etc..
- Collection of cells
- Radio-immuno tests
- Ligand binding test

Material:

- Glass fibre filter
- 100% borosilicated glass fibre
- Binder: free

Technical Specification	Values	Units	Methods
Weight	52	g/sqm	150536
Thickness	0.26	mm	
Filtration Velocity		sec.	DIN 53137
	100	sec.	Herzberg
Pore size	1.2	µm	
Dry bursting strength	2.2	psi	in-house

Disclaimer: VWR wants to state that these specifications can only be regarded as general information and will not discharge the user from his obligation to ensure the product is suitable for the intended use - If there's any doubt, ask us.

Date of issue: **September 2012**

VWR International bvba/sprl, Researchpark Haasrode 2020, Geldenaaksebaan 464, 3001 Leuven, Belgium

<http://www.vwr.com>



Version: 2.0