



Certificate of Conformity

VWR CALIBRATION WEIGHT SETS IN WOODEN BOXES, CLASS F2

VWR International certifies that the products specified conform to all specifications and quality criteria mentioned in this document.

VWR International certifie que les produits mentionnés respectent toutes les spécifications et critères de qualité mentionnés dans ce document.

VWR International bestätigt hiermit, daß die angegebenen Produkte alle in diesem Dokument erwähnten Spezifikationen und Qualitätskriterien erfüllen.

Catalog No.	Weight set	No. of weights	Total weight (g)	Pk
611-8662	1 mg - 500 mg	12	1.11	1
611-8663	1 mg - 50 g	20	111.11	1
611-8664	1 mg - 100 g	21	211.11	1
611-8665	1 mg - 200 g	23	611.11	1
611-8666	1 mg - 500 g	24	1111.11	1
611-8667	1 mg - 1 kg	25	2111.11	1
611-8668	1 mg - 2 kg	27	6111.11	1
611-8669	1 mg - 5 kg	28	11111.11	1
611-8670	1 mg - 10 kg	29	21111.11	1
611-8671	1 g - 50 g	8	110	1
611-8672	1 g - 100 g	9	210	1
611-8673	1 g - 200 g	11	610	1
611-8674	1 g - 500 g	12	1110	1
611-8675	1 g - 1 kg	13	2110	1
611-8676	1 g - 2 kg	15	6110	1
611-8677	1 g - 5 kg	16	11110	1
611-8678	1 g - 10 kg	17	21110	1

Other Characteristics:

- Wire weights up to 500 mg
- Identification number guarantees traceability for weights from 1 g upwards
- Available without COFRAC certificate
- Conform to OILM R111 recommendations.

Disclaimer: VWR wants to state that this declaration will not discharge the user from his obligation to ensure the product is suitable for the intended use. - If there is any doubt ask us.

Date of issue: **March 2015**