

PYREX[®] Laboratory Glassware

Fuel your expertise with the lab glassware trusted
by scientists for more than 100 years.





Table of Contents

Exceptional strength and stability. Brilliant transparency. Superior resistance to chemicals, contaminants, and drastic changes in temperature. These are among the key properties that enable PYREX® glass to endure even the toughest laboratory environments, making PYREX scientific glassware the indispensable partner for chemists and research scientists for more than 100 years.

PYREX® Glassware	3
PYREX VISTA™ Glassware	5
PYREXPLUS® Glassware	6
PYREX Low Actinic Glassware	7

PYREX[®] Glassware

Of the hundreds of commercial glasses produced, Code No. 7740 borosilicate glass comes closest to being the ideal glass for most laboratory applications. With proper care, it will withstand nearly all temperatures used in normal laboratory use. It is highly resistant to chemical attack. Its low coefficient of expansion allows it to be manufactured to suit a variety of wall thicknesses, tailored to different purposes. This mechanical strength is not solely dependent on wall thickness but is also influenced by the glass formula and type, ensuring durability while maintaining reliable heat resistance.



Ordering Information

Description	Volumes	VWR Cat. No.	Corning Cat. No.
Griffin Low Form Beaker	10 mL - 4L	13912-069	1000-10
		13912-091	1000-20
		13912-127	1000-30
		13912-149	1000-50
		13912-160	1000-100
		13912-182	1000-150
		13912-207	1000-250
		13912-229	1000-400
		13912-240	1000-600
		13912-262	1000-800
		13912-284	1000-1L
		13912-320	1000-2L
		89090-522	1000-3L
		89090-524	1000-4L
Round Media Storage Bottles, with Screw Cap	25 mL - 10L	16157-071	1395-25
		16157-068	1395-50
		16157-103	1395-100
		89427-656	1395-150
		16157-136	1395-250
		16157-169	1395-500
		89427-658	1395-750
		16157-191	1395-1L
		16157-227	1395-2L
		89427-655	1395-3X
		16157-250	1395-5L
		16157-282	1395-10L
Double Metric Scale Cylinders, Class A, To Deliver	25 mL - 2L	89090-496	3023-25
		89090-504	3023-50
		89090-488	3023-100
		89090-492	3023-250
		89090-500	3023-500
		22877-030	3023-1L
		22877-032	3023-2L

PYREX[®] Glassware (continued)

Ordering Information

Description	Volumes	VWR Cat. No.	Corning Cat. No.
Narrow Mouth Erlenmeyer Flasks with Heavy Duty Rim	10 mL - 6L	89091-116	4980-10
		29136-026	4980-25
		29136-037	4980-50
		29136-048	4980-125
		29136-060	4980-250
		89090-894	4980-300
		89091-420	4980-500
		29136-106	4980-1L
		470211-488	4980-2L
		89090-738	4980-4L
		89090-942	4980-6L
Class A Volumetric Flasks with PYREX Standard Taper Stopper	1 mL - 6L	89091-218	5640-1
		89091-300	5640-2
		89091-480	5640-5
		89091-102	5640-10
		89090-818	5640-25
		89090-912	5640-50
		89090-606	5640-100
		89090-680	5640-200
		89090-770	5640-250
		89090-752	5640-500
		89090-714	5640-1L
		89090-848	5640-2L
		26307-462	5640-3L
		26307-464	5640-4L
		26307-466	5640-5L
		26307-468	5640-6L
Culture Tube, Rimless	0.5 mL - 70 mL	470211-678	9820-6
		60820-080	9820-10
		60820-104	9820-12
		60820-126	9820-13
		60820-181	9820-16
		60820-206	9820-16X
		470211-668	9820-16XX
		470211-670	9820-18
		470211-672	9820-20
		60820-265	9820-22
		60820-320	9820-25
		60820-321	9820-25X
Culture Tube, Screw Caps	9 mL - 170 mL	470202-414	9825-13
		60827-045	9825-16
		470188-784	9825-16X
		60827-089	9825-16XX
		60827-125	9825-20
		60827-147	9825-20X
		89091-462	9825-22
		470211-682	9825-25
		60827-205	9825-25X
		60827-207	9825-38

PYREX[®] VISTA[™] Glassware

PYREX VISTA glassware is an economical option for the customer who is willing to forgo the premium benefits of PYREX products. Manufactured to Corning/PYREX standards and priced competitively against comparable products, PYREX VISTA glassware offers a full range of products from beakers to pipets and is easily recognized by its blue graduations and novel marking spot.



Ordering Information


Description	Volumes	VWR Cat. No.	Corning Cat. No.
Beakers, Standard Low Griffin	10 mL - 1L	89090-426	70000-10
		89090-432	70000-20
		89090-436	70000-30
		89090-440	70000-50
		89090-424	70000-100
		89090-428	70000-150
		470211-578	70000-250
		89090-438	70000-400
		89090-442	70000-600
		89090-430	70000-1L
Single Metric Scale Graduated Cylinder, Class A with Funnel Top, To Deliver	10 mL - 500 mL	89090-450	70024-10
		89090-458	70024-25
		89090-466	70024-50
		89090-446	70024-100
		89090-454	70024-250
		89090-462	70024-500
Narrow Mouth Erlenmeyer Flask	125 mL - 1L	89090-316	70980-125
		89090-348	70980-250
		89090-386	70980-500
		89090-302	70980-1L
Class A Volumetric Flasks with Glass Standard Taper Stopper	10 mL - 2L	89090-116	70640-10
		89090-354	70640-25
		89090-392	70640-50
		89090-304	70640-100
		89090-334	70640-200
		89090-342	70640-250
		89090-380	70640-500
		89090-298	70640-1L
		89090-330	70640-2L
Test Tube with Beaded Rim	9 mL - 70 mL	89090-422	70800-13
		89090-318	70800-15
		89090-340	70800-16
		89090-366	70800-18
		89090-374	70800-20
		89090-408	70800-25
		89090-416	70800-25X
Culture Tube, Screw Cap, Phenolic Cap	9 mL - 55 mL	89090-420	70825-13
		89090-312	70825-16
		89090-322	70825-16X
		89090-338	70825-16XX
		89090-362	70825-20
		89090-372	70825-20X
		89090-406	70825-25

PYREXPLUS® Glassware

PYREXPLUS glassware is coated with a tough, transparent plastic vinyl. The coating, which is applied to the outside of the vessel, helps prevent exterior surface abrasion. It also helps minimize the loss of contents and helps contain glass fragments if the glass vessel is broken.

Ordering Information

Description	Volumes	VWR Cat. No.	Corning Cat. No.
PVC-Coated Round Media Storage Bottles, with Screw Cap	100 mL - 2L	16157-057	61626-250
		16157-055	61626-500
		16157-053	61626-1L
		16157-051	61626-2L
Single Metric Scale Cylinders, Protective Coating, Graduated, To Deliver	100 mL - 2L	89090-578	63024-100
		89090-582	63024-250
		89090-590	63024-500
		89090-580	63024-1L
		89090-586	63024-2L
Heavy Wall Filtering Flask, Protective Coating, Tubulation, Graduated	250 mL - 4L	29410-988	65340-250
		89091-694	65340-500
		89091-676	65340-1L
		89090-584	65340-2L
		89090-588	65340-4L
Class A Volumetric Flask, Protective Coating, with Standard Taper Stopper	100 mL - 2L	29618-751	65640-100
		89091-682	65640-250
		89091-692	65640-500
		89091-674	65640-1L
		89091-688	65640-2L
Squibb Pear-shaped Separatory Funnels, Protective Coating, PTFE Product Standard Stopcock, Glass Standard Taper Stopper	125 mL - 2L	89091-670	66402-125
		89091-696	66402-500
		30356-249	66402-1L
		30356-252	66402-2L





PYREX[®] Low Actinic Glassware

PYREX Low Actinic Labware is made from Code. No. 7740 glass with copper replacing the sodium in the glass matrix at the surface. The resulting product is as durable as the base glass. PYREX Low Actinic Labware was originally developed for work in the vitamin field, but it has found other uses in applications with chemicals sensitive to light in the 3,000 to 5,000 Angstrom range. Low actinic glassware is designed to reduce light transmission, especially UV and visible wavelengths, protecting sensitive samples such as chemicals, reagents, and biological materials from degradation caused by light exposure. By minimizing light-induced reactions, the glassware ensure the stability and integrity of stored or processed substances, which is critical for precision in experiments and long-term storage. PYREX Low Actinic glassware is ideal for use in chemistry, biology, pharmaceuticals, and photochemistry experiments.

Ordering Information

Description	Volumes	VWR Cat. No.	Corning Cat. No.
Low Actinic Round Media Storage Bottles, Graduated, with Screw Cap	25 mL - 10L	33500-510	51395-25
		33500-512	51395-50
		6157-067	51395-100
		16157-065	51395-250
		16157-063	51395-500
		16157-061	51395-1L
		89090-098	51395-2L
		89090-104	51395-5L
		80080-920	51395-10L
Single Metric Scale Cylinders, Lifetime Red™, Graduated, To Contain	100 mL - 2L	89090-636	3042-10
		89090-624	3042-100
		89090-888	3042-500
		24710-522	3042-1L
		24710-544	3042-2L
Class A Low Actinic Volumetric Flask with PYREX Standard Taper Stopper	10 mL - 2L	89091-618	55640-10
		89090-540	55640-25
		89090-542	55640-50
		89090-530	55640-100
		89090-534	55640-200
		89090-538	55640-250
		29620-370	55640-500
		29620-372	55640-1L
		66025-530	55640-2L
Class A Lifetime Red Volumetric Flask with Standard Taper Stopper, To Contain	25 mL - 2L	89090-810	5660-25
		89090-902	5660-50
		89090-596	5660-100
		89090-670	5660-200
		89090-686	5660-250
		89090-742	5660-500
		89090-704	5660-1L
		89090-836	5660-2L



Warranty/Disclaimer: Unless otherwise specified, all products are for research use or general laboratory use only.* Not intended for use in diagnostic or therapeutic procedures. Not for use in humans. These products are not intended to mitigate the presence of microorganisms on surfaces or in the environment, where such organisms can be deleterious to humans or the environment. Corning Life Sciences makes no claims regarding the performance of these products for clinical or diagnostic applications. *For a listing of US medical devices, regulatory classifications or specific information on claims, visit www.corning.com/resources.

Corning's products are not specifically designed and tested for diagnostic testing. Many Corning products, though not specific for diagnostic testing, can be used in the workflow and preparation of the test at the customers discretion. Customers may use these products to support their claims. We cannot make any claims or statements that our products are approved for diagnostic testing either directly or indirectly. The customer is responsible for any testing, validation, and/or regulatory submissions that may be required to support the safety and efficacy of their intended application.

CORNING

For additional product or technical information, visit vwr.com/corning, call 1.800.932.5000, or contact your Avantor representative.



vwr.com

Any prices, product, and/or services details are current when published and subject to change. Offers subject to product availability. Void where prohibited by law or company policy. Any claims concerning products manufactured by third parties are the sole claims of the manufacturer and not those of Avantor or its affiliates. All products are subject to the terms and conditions of sale of Avantor or its affiliates, including the applicable limitations of liability contained therein. The information herein should not be considered a warranty of any kind. Please carefully review any product information sheets and safety information for the product. Trademarks are owned by Avantor Inc. or its affiliates, unless otherwise noted. | Visit vwr.com to view our privacy policy and additional disclaimers. © 2026 Avantor Inc. All rights reserved.