

Product Information Data Sheet

VWR Catalogue Number	See List
Country of Origin	Manufactured in Czech Republic
Date of Issue (yyyy-mm-dd)	2026-07-21

Quality System Compliance

This certificate confirms that the product has been sourced from a supplier or manufacturer who has declared their Quality Management System (QMS) complies with **ISO 9001 & ISO 14001**.

QC Testing

Representative products samples are collected and inspected in accordance with current applicable product specifications.

Product Specifications

Material	Cardboard CAMA white/grey, 700 gsm, 1.0 mm, FSC recycled 100% + Arctic paper G print 130 gsm.
Sterilization	Not Applicable.
ATP Assay	Not applicable.
DNase & RNase Free	Not Applicable.
BSE/TSE	Not applicable.
Cytotoxicity	Not applicable.
Non-Pyrogenic Statement	Not applicable.
Latex Statement	Not applicable.
BPA Statement	Not applicable.
DEHP Statement	Not applicable.
RoHS	Not applicable.

REACH Statement	Not applicable.
Storage Conditions	Store in Humidity max 80%.
Shelf Life	Not relevant.

Product Family

Mechanical cryogenic freezer boxes with or without dividers

VWR Catalog Number	Description
479-2220	Cryobox without divider, waterproof, White, 127×127×73.0 mm
479-2215	Cryobox without divider, waterproof, White, 133×133×50 mm
479-2217	Cryobox with divider (height 30 mm), waterproof, with holes in the bottom, 25 (5×5) Compartments, White, 75×75×50 mm
479-2213	Cryobox with divider (height 30 mm), waterproof, with holes in the bottom, 81 (9×9) Compartments, White, 133×133×50 mm
479-2227	Cryobox with divider (height 20 mm), water repellent, for 0.2 ml PCR tubes, 81 (9×9) Compartments, White, 76×76×24.4 mm
479-2228	Cryobox with divider (height 20 mm), water repellent, for 0.2 ml PCR tubes, 196 (14×14) Compartments, White, 133×133×24.4 mm
479-0249	81 Compartments, 127×127×47.6 mm

Fibreboard dividers

VWR Catalog Number	Description
479-2225	16 Compartments, 125×125×25.4 mm
479-2224	25 Compartments, 133×133× 25.4 mm
479-2223	49 Compartments, 125×125×25.4 mm
479-2219	64 Compartments, 133×133× 25.4 mm
479-2221	81 Compartments, 125×125×25.4 mm
479-2216	100 Compartments, 133×133×25.4 mm

Disclaimer: VWR'S Quality Management System (QMS) is ISO 9001:2015 certified as a wholesale distributor of finished products, ensuring effective oversight of finished product handling, traceability and distribution processes. VWR is not a manufacturer and does not engage in manufacturing, co- manufacturing, or processing activities such as design, packaging, labelling, testing, or quality release. Consequently, the supplier or manufacturer's QMS certification may differ from VWR's.

VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. This product is intended for research use only. It is not for use in diagnostic procedures, therapeutic applications, or any human or clinical use. This document has been produced electronically and is valid without a signature.

Version History

Sl.no	Version number	Description	Month & Year
1	1.0.0	Created PIDS	07-2026