

Certificate of Conformity

VWR Catalogue Number	VWRA548-2672
Description	VWR® Snap Ring Vials with write-on spot, ND11, 2 mL capacity, 11.6 mm diameter, 32 mm height, without cap, Amber, snap cap type
Lot Number	20240429
Date of Manufacture (yyyy-mm-dd)	2024-04-29
Country of Origin	Manufactured in China
Date of Issue (yyyy-mm-dd)	2025-05-26

Quality System Compliance

Products are manufactured under the **ISO 9001:2015 & GB/T 19001-2016** standard. Products are inspected and controlled through whole production processing in accordance with current applicable product specifications and QC SOP.

QC Testing

Representative products samples are collected and inspected in accordance with current applicable QC requirements.

Test/Procedure

Inside Surface Hydrolytic Resistance
Content of B₂O₃
Coefficient of Thermal Expansion
Hydrolytic resistance of glass grains at 121 °C
Thermal shock endurance resistance
Stress
Leaching amount of As、Sb、Pb、Cd

Result

Accord with HC1
Not be less than 8%(g/g), Accord with USP I
(3.5~6.1)*10⁻⁶ K⁻¹ (20°C~300°C), Accord with USP I.
Accord with level 1 by USP I.
Pass
Pass
As≤0.2mg/L; Sb≤0.7mg/L; Pb≤1.0mg/L;
Cd≤0.25mg/L.

Product Specifications

Material Tube: 5.0 Borosilicate glass (Type I)

Sterilization	Not Applicable.
ATP Assay	Not Applicable.
DNase & RNase Free	Not Applicable.
BSE/TSE	No use of any raw material produced from or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE) according to EMEA 410/01 revision 3.
Cytotoxicity	Not Applicable.
Non-Pyrogenic Statement	Not Applicable.
Latex Statement	This product is not made from natural rubber latex.
BPA Statement	Not Applicable.
DEHP Statement	Not Applicable.
RoHS	Not Applicable.
REACH Statement	Not Applicable.
Storage Conditions	Store at normal temperature.
Shelf Life	Not Applicable.

Disclaimer: VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. The purpose of the product is for use in laboratory only.
This document has been produced electronically and is valid without a signature.