

SAFETY DATA SHEET

according to Regulation (EC) No. 1907/2006

Revision Date 05.06.2017

Version 1.4

SECTION 1. Identification of the substance/mixture and of the company/undertaking

1.1 Product identifier

Catalogue No.	199002
Product name	Buffer solution (potassium dihydrogen phosphate/di-sodium hydrogen phosphate), traceable to SRM from NIST and PTB pH 7.00 (25°C) CertiPUR®

REACH Registration Number This product is a mixture. REACH Registration Number see section 3.

1.2 Relevant identified uses of the substance or mixture and uses advised against

Identified uses	Reagent for development and research For additional information on uses please refer to the Merck Chemicals portal (www.merckgroup.com).
-----------------	---

1.3 Details of the supplier of the safety data sheet

Responsible Department	LS-QHC * e-mail: prodsafe@merckgroup.com
Regional representation	Merck Chemicals Ltd * Boulevard Industrial Park * Padge Road * Beeston * Nottingham * NG9 2JR * Tel. 01159 430840 * information@merckchem.co.uk .

1.4 Emergency telephone number	+49 (0) 6151 722440
--------------------------------	---------------------

SECTION 2. Hazards identification

2.1 Classification of the substance or mixture

This mixture is not classified as dangerous according to European Union legislation.

SAFETY DATA SHEET

according to Regulation (EC) No. 1907/2006

Catalogue No.	199002
Product name	Buffer solution (potassium dihydrogen phosphate/di-sodium hydrogen phosphate), traceable to SRM from NIST and PTB pH 7.00 (25° C) CertiPUR®

2.2 Label elements

Labelling (REGULATION (EC) No 1272/2008)

Not a hazardous substance or mixture according to Regulation (EC) No. 1272/2008.

2.3 Other hazards

None known.

SECTION 3. Composition/information on ingredients

Chemical nature	Aqueous solution of inorganic compounds.
-----------------	--

3.1 Substance

Not applicable

3.2 Mixture

Remarks	No disclosure requirement according to Regulation (EC) No. 1907/2006.
---------	---

SECTION 4. First aid measures

4.1 Description of first aid measures

After inhalation: fresh air.

In case of skin contact: Take off immediately all contaminated clothing. Rinse skin with water/shower.

After eye contact: rinse out with plenty of water. Remove contact lenses.

After swallowing: make victim drink water (two glasses at most). Consult doctor if feeling unwell.

4.2 Most important symptoms and effects, both acute and delayed

We have no description of any toxic symptoms.

SAFETY DATA SHEET

according to Regulation (EC) No. 1907/2006

Catalogue No.	199002
Product name	Buffer solution (potassium dihydrogen phosphate/di-sodium hydrogen phosphate), traceable to SRM from NIST and PTB pH 7.00 (25° C) CertiPUR®

4.3 Indication of any immediate medical attention and special treatment needed

No information available.

SECTION 5. Firefighting measures

5.1 Extinguishing media

Suitable extinguishing media

Use extinguishing measures that are appropriate to local circumstances and the surrounding environment.

Unsuitable extinguishing media

For this substance/mixture no limitations of extinguishing agents are given.

5.2 Special hazards arising from the substance or mixture

Not combustible.

Ambient fire may liberate hazardous vapours.

5.3 Advice for firefighters

Special protective equipment for firefighters

In the event of fire, wear self-contained breathing apparatus.

SECTION 6. Accidental release measures

6.1 Personal precautions, protective equipment and emergency procedures

Advice for non-emergency personnel: Do not breathe vapours, aerosols. Evacuate the danger area, observe emergency procedures, consult an expert.

Advice for emergency responders:

Protective equipment see section 8.

6.2 Environmental precautions

No special precautionary measures necessary.

6.3 Methods and materials for containment and cleaning up

SAFETY DATA SHEET

according to Regulation (EC) No. 1907/2006

Catalogue No.	199002
Product name	Buffer solution (potassium dihydrogen phosphate/di-sodium hydrogen phosphate), traceable to SRM from NIST and PTB pH 7.00 (25° C) CertiPUR®

Observe possible material restrictions (see sections 7 and 10).

Take up with liquid-absorbent material (e.g. Chemizorb®). Dispose of properly. Clean up affected area.

6.4 Reference to other sections

Indications about waste treatment see section 13.

SECTION 7. Handling and storage

7.1 Precautions for safe handling

Advice on safe handling

Observe label precautions.

Hygiene measures

Change contaminated clothing. Wash hands after working with substance.

7.2 Conditions for safe storage, including any incompatibilities

Storage conditions

Tightly closed.

Recommended storage temperature see product label.

7.3 Specific end use(s)

Apart from the uses mentioned in section 1.2 no other specific uses are stipulated.

SECTION 8. Exposure controls/personal protection

8.1 Control parameters

Contains no substances with occupational exposure limit values.

8.2 Exposure controls

Engineering measures

SAFETY DATA SHEET

according to Regulation (EC) No. 1907/2006

Catalogue No.	199002
Product name	Buffer solution (potassium dihydrogen phosphate/di-sodium hydrogen phosphate), traceable to SRM from NIST and PTB pH 7.00 (25° C) CertiPUR®

Technical measures and appropriate working operations should be given priority over the use of personal protective equipment.

See section 7.1.

Individual protection measures

Protective clothing needs to be selected specifically for the workplace, depending on concentrations and quantities of the hazardous substances handled. The chemical resistance of the protective equipment should be enquired at the respective supplier.

Eye/face protection

Safety glasses

Hand protection

not required

Respiratory protection

Not required; except in case of aerosol formation.

Environmental exposure controls

No special precautionary measures necessary.

SECTION 9. Physical and chemical properties

9.1 Information on basic physical and chemical properties

Form	liquid
Colour	colourless
Odour	odourless
Odour Threshold	Not applicable

SAFETY DATA SHEET

according to Regulation (EC) No. 1907/2006

Catalogue No.	199002
Product name	Buffer solution (potassium dihydrogen phosphate/di-sodium hydrogen phosphate), traceable to SRM from NIST and PTB pH 7.00 (25° C) CertiPUR®

pH	7.0 at 25 °C
----	-----------------

Melting point	-5 °C
---------------	-------

Boiling point/boiling range	109 °C at 1,013 hPa
-----------------------------	------------------------

Flash point	Not applicable
-------------	----------------

Evaporation rate	No information available.
------------------	---------------------------

Flammability (solid, gas)	No information available.
---------------------------	---------------------------

Lower explosion limit	No information available.
-----------------------	---------------------------

Upper explosion limit	No information available.
-----------------------	---------------------------

Vapour pressure	No information available.
-----------------	---------------------------

Relative vapour density	No information available.
-------------------------	---------------------------

Density	1.01 g/cm ³ at 20 °C
---------	------------------------------------

Relative density	No information available.
------------------	---------------------------

Water solubility	at 20 °C soluble
------------------	---------------------

Partition coefficient: n-octanol/water	No information available.
--	---------------------------

Auto-ignition temperature	No information available.
---------------------------	---------------------------

SAFETY DATA SHEET

according to Regulation (EC) No. 1907/2006

Catalogue No.	199002
Product name	Buffer solution (potassium dihydrogen phosphate/di-sodium hydrogen phosphate), traceable to SRM from NIST and PTB pH 7.00 (25° C) CertiPUR®

Decomposition temperature	No information available.
Viscosity, dynamic	No information available.
Explosive properties	Not classified as explosive.
Oxidizing properties	none

9.2 Other data

none

SECTION 10. Stability and reactivity

10.1 Reactivity

See section 10.3

10.2 Chemical stability

The product is chemically stable under standard ambient conditions (room temperature) .

10.3 Possibility of hazardous reactions

Violent reactions possible with:

The generally known reaction partners of water.

10.4 Conditions to avoid

no information available

10.5 Incompatible materials

no information available

10.6 Hazardous decomposition products

no information available

SAFETY DATA SHEET

according to Regulation (EC) No. 1907/2006

Catalogue No.	199002
Product name	Buffer solution (potassium dihydrogen phosphate/di-sodium hydrogen phosphate), traceable to SRM from NIST and PTB pH 7.00 (25° C) CertiPUR®

SECTION 11. Toxicological information

11.1 Information on toxicological effects

Mixture

Acute oral toxicity

This information is not available.

Acute inhalation toxicity

This information is not available.

Acute dermal toxicity

This information is not available.

Skin irritation

This information is not available.

Eye irritation

This information is not available.

Sensitisation

This information is not available.

Germ cell mutagenicity

This information is not available.

Carcinogenicity

This information is not available.

Reproductive toxicity

This information is not available.

Teratogenicity

This information is not available.

Specific target organ toxicity - single exposure

This information is not available.

Specific target organ toxicity - repeated exposure

This information is not available.

SAFETY DATA SHEET

according to Regulation (EC) No. 1907/2006

Catalogue No.	199002
Product name	Buffer solution (potassium dihydrogen phosphate/di-sodium hydrogen phosphate), traceable to SRM from NIST and PTB pH 7.00 (25° C) CertiPUR®

Aspiration hazard

This information is not available.

11.2 Further information

Hazardous properties cannot be excluded but are unlikely when the product is handled appropriately.

Handle in accordance with good industrial hygiene and safety practice.

SECTION 12. Ecological information

Mixture

12.1 Toxicity

No information available.

12.2 Persistence and degradability

No information available.

12.3 Bioaccumulative potential

No information available.

12.4 Mobility in soil

No information available.

12.5 Results of PBT and vPvB assessment

Substance(s) in the mixture do(es) not meet the criteria for PBT or vPvB according to Regulation (EC) No 1907/2006, Annex XIII, or a PBT/vPvB assessment was not conducted.

12.6 Other adverse effects

Discharge into the environment must be avoided.

SAFETY DATA SHEET

according to Regulation (EC) No. 1907/2006

Catalogue No.	199002
Product name	Buffer solution (potassium dihydrogen phosphate/di-sodium hydrogen phosphate), traceable to SRM from NIST and PTB pH 7.00 (25° C) CertiPUR®

SECTION 13. Disposal considerations

Waste treatment methods

Waste material must be disposed of in accordance with the national and local regulations. Leave chemicals in original containers. No mixing with other waste. Handle uncleaned containers like the product itself.

See www.retrologistik.com for processes regarding the return of chemicals and containers, or contact us there if you have further questions.

Notice Directive on waste 2008/98/EC.

SECTION 14. Transport information

Land transport (ADR/RID)

14.1 - 14.6	Not classified as dangerous in the meaning of transport regulations.
-------------	--

Inland waterway transport (ADN)

Not relevant

Air transport (IATA)

14.1 - 14.6	Not classified as dangerous in the meaning of transport regulations.
-------------	--

Sea transport (IMDG)

14.1 - 14.6	Not classified as dangerous in the meaning of transport regulations.
-------------	--

14.7 Transport in bulk according to Annex II of MARPOL 73/78 and the IBC Code

Not relevant

SAFETY DATA SHEET

according to Regulation (EC) No. 1907/2006

Catalogue No.	199002
Product name	Buffer solution (potassium dihydrogen phosphate/di-sodium hydrogen phosphate), traceable to SRM from NIST and PTB pH 7.00 (25° C) CertiPUR®

SECTION 15. Regulatory information

15.1 Safety, health and environmental regulations/legislation specific for the substance or mixture

EU regulations

Major Accident Hazard	SEVESO III
Legislation	Not applicable

Regulation (EC) No 1005/2009 on substances that deplete the ozone layer	not regulated
---	---------------

Regulation (EC) No 850/2004 of the European Parliament and of the Council of 29 April 2004 on persistent organic pollutants and amending Directive 79/117/EEC	not regulated
---	---------------

Substances of very high concern (SVHC)	This product does not contain substances of very high concern according to Regulation (EC) No 1907/2006 (REACH), Article 57 above the respective regulatory concentration limit of ≥ 0.1 % (w/w).
--	--

National legislation

Storage class	10 - 13
---------------	---------

15.2 Chemical safety assessment

For this product a chemical safety assessment was not carried out.

SECTION 16. Other information

Training advice

Provide adequate information, instruction and training for operators.

SAFETY DATA SHEET

according to Regulation (EC) No. 1907/2006

Catalogue No. 199002

Product name Buffer solution (potassium dihydrogen phosphate/di-sodium hydrogen phosphate), traceable to SRM from NIST and PTB pH 7.00 (25° C)
CertiPUR®

Key or legend to abbreviations and acronyms used in the safety data sheet

Used abbreviations and acronyms can be looked up at www.wikipedia.org.

The information contained herein is based on the present state of our knowledge. It characterises the product with regard to the appropriate safety precautions. It does not represent a guarantee of any properties of the product.