

PDS No. 189170	PRODUCT DATA SHEET	Page 1 of 1
Revision 02	Tube, 7 ml, PS, Screw Cap	 greiner bio-one
	Greiner Item-No. 189170	

1.	Description / Specification	
1.1	Description	PS Tube with flat bottom, 7 ml, Screw Cap
1.2	Dimensions	<u>Tube</u> : 48,6 x 17,6 mm <u>Tube with cap</u> : 50,2 x 22,6 mm
1.3	Volume	7 ml
1.4	Material / Resin	<u>Tube</u> : PS (Polystyrene), free of heavy metal <u>Cap</u> : HDPE (High Density Polyethylene), free of heavy metal
1.5	Colour	<u>Tube</u> : clear <u>Cap</u> : white
1.6	Sterilization	No
1.7	Quality Control	- <u>Raw Material-Control</u> : physical testing - <u>Product-Control</u> : testing of attributive and variable characteristics in accordance with the valid specification
1.8	Other Information	For single use only

2.	Features	
2.1	Basic features	-
2.2	Temperature range	-20°C to +60°C
2.3	Autoclavability	No
2.4	Centrifugation	7200 x g: fixed-angle rotor
2.5	Chemical Resistance	See homepage: www.gbo.com/bioscience →Products →Literature →Technical Information→Chemical Resistance of Resins
2.6	Shelf life	N/A
2.7	Other Information	-

3.	Packaging	
3.1	Pieces / Bag	700
3.2	Pieces / Box	700
3.3	Lot-No.	YYWWXX (year, week, day, shift)
3.4	Other Information	-

4.	Other Information	

Data Sheet subject to change without notice!

Prior Issue	Drawn	Approved	Released	CONFIDENTIAL : Information contained in this document or drawing is confidential and proprietary to Greiner Bio-One GmbH. This document may not be reproduced for any reason without written permission from Greiner Bio-One GmbH. All rights of design, invention, and copyright are reserved.
Revision 01	Date 16 December 2009	Date 17 December 2009	Date 17 December 2009	
Date 22.05.2007	Name S. Kaelberer	Name Dr. L. Marchetti	Name A. Schulz	