

Certificate of Conformity

VWR Catalogue Number	See List
Description	VWR®, Isolation Gowns
Country of Origin	Manufactured in China
Date of Issue (yyyy-mm-dd)	2024-03-18

45 g/m²

Catalogue Number	Colour	W x L [mm]
113-8316	Blue	1150 × 1360
113-8317	Green	1150 × 1360
113-8319	Pink	1150 × 1360
113-8318	Red	1150 × 1360
113-8315	White	1150 × 1360
113-8320	Yellow	1150 × 1360

60 g/m²

Catalogue Number	Colour	W x L [mm]
113-8322	Blue	1150 × 1360
113-8323	Green	1150 × 1360
113-8325	Pink	1150 × 1360
113-8324	Red	1150 × 1360
113-8321	White	1150 × 1360
113-8326	Yellow	1150 × 1360

Quality System Compliance

Products are manufactured under the **ISO 13485:2016** standard. Products are inspected and controlled through whole production processing in accordance with current applicable product specifications and QC SOP.

QC Testing

Representative products samples are collected and inspected in accordance with current applicable QC requirements. Inspected according to AQL level II 4.0 sampling plan per MIL-STD-105E.

45 g/m²

	Test Method	MD	CD
Break Tensile Properties [N]	ISO13934-1	> 50	> 40
Elongation rate [%]	ISO13934-1	> 30	> 30

60 g/m²

	Test Method	MD	CD
Break Tensile Properties [N]	ISO13934-1	> 70	> 55
Elongation rate [%]	ISO13934-1	> 30	> 30

Product Specifications

Material	Non-woven Spunbond Meltblown Spunbond.
Sterilization	Non-Sterile.
ATP Assay	Not Applicable.
DNase & RNase Free	Not Applicable.
BSE/TSE	No use of any raw material produced from, or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE).
Cytotoxicity	Not Applicable.
Non-Pyrogenic Statement	Not Applicable.
Latex Statement	This product is not made from natural rubber latex.

BPA Statement	Not Applicable.
DEHP Statement	Not Applicable.
RoHS	Not Applicable.
REACH Statement	Not Applicable.
Storage Conditions	Not Applicable.
Shelf Life	Not Applicable.

Disclaimer: VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. The purpose of the product is for use in laboratory only.
This document has been produced electronically and is valid without a signature.

Version History

Sl.no	Version number	Description	Month & Year
1	1.0	Created General CoQ using new layout	03-2024