

Certificate of Conformity

VWR Catalogue Number	548-2644A
Description	VWR® Short Thread Vials ND9 with integrated micro-insert, 32.0 × 11.6 mm, 0.2 mL capacity, clear
Country of Origin	Manufactured in Germany
Date of Issue (yyyy-mm-dd)	2024-05-20

Quality System Compliance

Products are manufactured under the **ISO 9001:2015** standard. Products are inspected and controlled through whole production processing in accordance with current applicable product specifications and QC SOP.

QC Testing

Representative products samples are collected and inspected in accordance with current applicable product specifications.

Visual Inspection: No deformation, contamination, enclosure.

Product Specifications

Material	Clear glass, 1 st hydrolytic class.
Sterilization	Non-Sterile.
ATP Assay	The product does not contain ATP.
DNase & RNase Free	This product is free of any detectable DNase/RNase contamination.
BSE/TSE	No use of any raw material produced from, or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE).
Cytotoxicity	This product is free of cytotoxicity effects.
Non-Pyrogenic Statement	Not applicable.

Latex Statement	This product is not made from natural rubber latex.
BPA Statement	Bisphenols are not used in the manufacture of the raw material and are not expected to be present.
DEHP Statement	Please see RoHS information below.
RoHS	Product does not contain hazardous substance restricted in Annex II of EU RoHS directive 2015/863/EU.
REACH Statement	To the best of our knowledge material used for the production does not contain more than 0.1 % (w/w) of any substance of very high concern.
Storage Conditions	Dry, protected from light and standard warehouse temperature.
Shelf Life	2 years, conditions: store clean, dry, and protected from light at normal storage temperature.

Disclaimer: VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. The purpose of the product is for use in laboratory only.
This document has been produced electronically and is valid without a signature.

Version History

Sl.no	Version number	Description	Month & Year
1	1.0	Created General CoC using new layout	05-2024