

# Certificate of Conformity

|                                   |                        |
|-----------------------------------|------------------------|
| <b>VWR Catalogue Number</b>       | See List               |
| <b>Description</b>                | VWR Macro Pipette Tips |
| <b>Country of Origin</b>          | Manufactured in USA    |
| <b>Date of Issue (yyyy-mm-dd)</b> | 2023-06-30             |

| Catalogue Number | Description                                                      |
|------------------|------------------------------------------------------------------|
| 613-7096         | Pipette Tips, 5 mL Eppendorf Fit, Bulk, Non-sterile              |
| 732-3914         | Filter Tips, 5 mL Eppendorf Fit, Bulk, Non-sterile               |
| 732-3915         | Filter Tips, 5 mL Eppendorf Fit, Racked, Sterile                 |
| 613-7097         | Pipette Tips, 5 mL Eppendorf Fit, Racked, Sterile                |
| 613-7094         | Pipette Tips, 5 mL Gilson Fit, Bulk, Non-sterile                 |
| 732-3912         | Filter Tips, 5 mL Gilson Fit, Bulk, Non-sterile                  |
| 732-3913         | Filter Tips, 5 mL Gilson Fit, Racked, Sterile                    |
| 613-7095         | Pipette Tips, 5 mL Gilson Fit, Racked, Sterile                   |
| 613-7100         | Pipette Tips, 10 mL Gilson/ Eppendorf Fit, Bulk, Non-sterile     |
| 732-3916         | Filter Tips, 10 mL Gilson/ Eppendorf Fit, Bulk, Non-sterile      |
| 732-3917         | Filter Tips, 10 mL Gilson/ Eppendorf Fit, Racked, Sterile        |
| 613-7101         | Pipette Tips, 10 mL Gilson/ Eppendorf Fit, Racked, Sterile       |
| 613-7098         | Pipette Tips, 10 mL Thermo Finn / Socorex Fit, Bulk, Non-sterile |
| 613-7099         | Pipette Tips, 10 mL Thermo Finn / Socorex Fit, Racked, Sterile   |

## Quality System Compliance

Products are manufactured under the **ISO 9001:2015** standard. Products are inspected and controlled through whole production processing in accordance with current applicable product specifications and QC SOP.

## QC Testing

Representative products samples are collected and inspected in accordance with current applicable QC requirements.

VWR International bv, Haasrode Researchpark 2020, Geldenaaksebaan 464, B-3001 Leuven, Belgium - <http://www.vwr.com>

## Product Specifications

The product is in conformance with the requirements of the following standards and regulations:

|                                |                                                                                                                                                                                                                                                                                                                                                                       |
|--------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Material</b>                | Virgin polypropylene.                                                                                                                                                                                                                                                                                                                                                 |
| <b>Sterilization</b>           | Products are treated with Electron-beam (E-beam) irradiation with a specified dose range of 12-25 kGy (Which is not applicable for non-Sterile products).                                                                                                                                                                                                             |
| <b>ATP Assay</b>               | The product is ATP free.                                                                                                                                                                                                                                                                                                                                              |
| <b>DNase &amp; RNase Free</b>  | The above products are free of any detectable DNase/RNase contamination.                                                                                                                                                                                                                                                                                              |
| <b>DNA Free</b>                | Product is DNA free.                                                                                                                                                                                                                                                                                                                                                  |
| <b>BSE/TSE</b>                 | No use of any raw material produced from or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE).           |
| <b>Cytotoxicity</b>            | Testing has been conducted and have been shown to be non-toxic.                                                                                                                                                                                                                                                                                                       |
| <b>Non-Pyrogenic Statement</b> | The acceptance level for product is 0.06 EU/mL (LAL Gel Clot Method).                                                                                                                                                                                                                                                                                                 |
| <b>Latex Statement</b>         | The above products are not made from natural rubber latex.                                                                                                                                                                                                                                                                                                            |
| <b>BPA Statement</b>           | Bisphenols are not used in the manufacture of the raw material and are not expected to be present.                                                                                                                                                                                                                                                                    |
| <b>DEHP Statement</b>          | Please see RoHS information below.                                                                                                                                                                                                                                                                                                                                    |
| <b>RoHS</b>                    | No substances (Lead, Cadmium, Mercury, Hexavalent Chromium (Cr6+), Poly Brominated Biphenyls (PBB), Poly Brominated Diphenyl ethers (PBDE), Bis(2-Ethyhexyl) Phthalate (DEHP), Benzyl Butyl Phthalate (BBP), Dibutyl Phthalate (DBP), Diisobutyl Phthalate (DIBP)) are used in manufacturing the raw materials and final product. No routinely analysis is performed. |

|                           |                                                                                                                                                                                                                  |
|---------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>REACH Statement</b>    | No substances as listed on >Substance of Very High Concern< candidate list published 2023-01-17 are used in manufacturing of the raw materials used to produce this product. No routinely analysis is performed. |
| <b>Storage Conditions</b> | Ambient Temperature.                                                                                                                                                                                             |
| <b>Shelf Life</b>         | 60 Months.                                                                                                                                                                                                       |

Disclaimer: VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. The purpose of the product is for use in laboratory only.  
This document has been produced electronically and is valid without a signature

## Version History

| Sl.no | Version number | Description                          | Month & Year |
|-------|----------------|--------------------------------------|--------------|
| 1     | 1.0            | Created General CoC using new layout | 06-2023      |
|       |                |                                      |              |