

Certificate of Conformity

VWR Catalogue Number	129-1552
Description	Autoclavable Biohazard Red Colour Bags with 50 µm thickness, 27 L Capacity, and 610×483 mm dimension
Lot Number	000964060300
Date of Manufacture (yyyy-mm-dd)	2022-04-18
Country of Origin	Manufactured in USA
Date of Issue (yyyy-mm-dd)	2022-08-25

Quality System Compliance

Products are manufactured under the **ISO 9001:2015** standard. Products are inspected and controlled through whole production processing in accordance with current applicable product specifications and QC SOP.

QC Testing

Representative products samples are collected and inspected in accordance with current applicable QC requirements.

Product Specifications

Material	Polypropylene.
Sterilization	Not Applicable. Products can be autoclaved at operating temperatures of 275 °F (135 °C).
ATP Assay	Not Applicable.
DNase & RNase Free	Not Applicable.
BSE/TSE	Resins and films used to manufacture this product may contain byproducts of bovine sourced raw material in their manufacture. These animal byproduct additives are processed under highly refined and controlled industrial chemical processes and at very high temperatures prior to purification and subsequent use as a polymer additive. The animal derived byproducts are not added during the process of manufacturing but originate either in the manufacture of the resins or in the additives used in the film formulations.

Cytotoxicity	Not Applicable.
Non-Pyrogenic Statement	Not Applicable.
Latex Statement	This product is not made from natural rubber latex.
BPA Statement	Not Applicable.
DEHP Statement	Please see RoHS information below.
RoHS	No substances (Lead, Cadmium, Mercury, Hexavalent Chromium (Cr6+), Poly Brominated Biphenyls (PBB), Poly Brominated Diphenyl ethers (PBDE), Bis(2-Ethyhexyl) Phthalate (DEHP), Benzyl Butyl Phthalate (BBP), Dibutyl Phthalate (DBP), Diisobutyl Phthalate (DIBP)) are used in manufacturing the raw materials and final product. No routinely analyse is performed.
REACH Statement	No substances as listed on >Substance of Very High Concern< candidate list published 2022-01-17 are used in manufacturing of the raw materials used to produce this product. No routinely analysis is performed.
Storage Conditions	Room temperature.

Disclaimer: VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. The purpose of the product is for use in laboratory only.
This document has been produced electronically and is valid without a signature.