



Product Specification Sheet

VWR CARBOYS, HIGH DENSITY POLYETHYLENE, WITH SCREW CAP, WITH SPIGOT

VWR International certifies that the specified products meet all VWR specifications and quality criteria for release. The products comply in every respect with the relevant particular specifications, drawings, relevant standards and regulations in force.

VWR International certifie que les produits mentionnés respectent toutes les spécifications et critères de qualité VWR. Les produits satisfont aux spécifications particulières, aux plans, aux normes et règlements en vigueur s'y rapportant.

VWR International bestätigt hiermit, daß die angegebenen Produkte alle VWR Spezifikationen und Qualitätskriterien zur Freigabe erfüllen. Die Produkte entsprechen in jeder Hinsicht den maßgeblichen Spezifikationen, den Zeichnungen sowie den relevanten Normen und gültigen Vorschriften.

Cat. No.	Capacity (L)	Fill Capacity (L)	Cap	Height With Cap(mm)	Length (mm)	Width (mm)	Spigot Extend (mm)	Neck Ø int (mm)	Neck Ø ext (mm)
216-0615	5	6.8	80	414	220	133	82	70.4	80.0
216-0622	10	13.0	80	490	274	164	52	70.4	80.0
216-0629	20	26.0	80	614	352	215	39	70.4	80
216-0635	20	26.0	120	640	352	215	39	110.7	120.0
216-0642	40	51.0	120	766	435	252	40	110.7	120.0

Bottle Material: High Density Polyethylene (HDPE)

Screw Cap Material: Polypropylene (PP)

Gasket Material: Platinum cured silicone

Spigot: Polypropylene (PP), DEHP free Flexible Polyvinyl Chloride (fPVC)

Properties:

- White
- Translucent
- Leakproof
- Scratch resistant
- ETO Stable
- Graduations to a +/- 5% accuracy
- Can be used with various magnetic mixers
- **Pressure and vacuum not recommended**
- **Do not autoclave!**

(*) observe chemical and temperature resistance



Product Specification Sheet

VWR CARBOYS, HIGH DENSITY POLYETHYLENE, WITH SCREW CAP, WITH SPIGOT

Temperature Resistance: -50°C to 120°C

- Do not maintain at extreme temperatures
- Do not expose to open flame, or heaters!

Chemical Resistance:

- Excellent chemical resistance at room temperature to most mild Acids, Bases/Alkali, and Aliphatic Alcohols.
- Good chemical resistance to most strong acids, aldehydes and their derivatives, esters, and aliphatic hydrocarbons.
- Very limited resistance to Ether, Aromatic and Halogenated Hydrocarbons, and ketones and derivatives.
- Do not use in combination with oxidizing acids and agents!
- Do not use for long-term storage of strong organic compounds
- Reduced chemical resistance at > 50°C.
- Comprehensive expected chemical resistance chart available upon request.
- Always test to specific chemical applications before use.

Product Compliance:

The above mentioned products, comply by composition with:

- The US FDA Title 21, CFR 177.1520 (c) 3.2a.
- US Pharmacopeia Class VI requirements
- The EU Directive 76/769/EC on the restriction on marketing and use of certain dangerous substances and preparations and following adaptations to technical progress, including EU directive 2005/84/EC on restrictions of phthalates.
- The European Directive 2002/79/EC and following amendments (2004/1/EC, 2004/19/EC, 2005/79/EC, 2007/19/EC, 2008/39/EC) on plastic materials intended to come into contact with foodstuffs
- The EU directive on Packaging and Packaging waste 94/62/EC-art 11

Disclaimer: VWR wants to state that these specifications can only be regarded as general information and will not discharge the user from his obligation to ensure the product is suitable for the intended use- If there is any doubt ask us.