



Product Information

1. Article number: 4098102
2. Product name: Discofix®
3. Product description: Multidirectional stopcocks with tubing for infusion therapy
4. List of materials: see attachments
5. Form of delivery: ☒ sterile ☐ non-sterile ☐ non-sterile/bulk
6. Information about sterilization: see attachments
7. Latex (natural rubber): ☐ yes ☒ no
8. DEHP (according to REACH regulation) ☐ yes ☒ no
9. Shelf life: ☐ 5 years ☒ 3 years ☐ not applicable
10. Instructions for use: ☐ yes ☒ not applicable
11. Bioburden: our non-sterile products do not exceed the Bioburden value of 200 CFU (colony forming unit)
12. Is the product CE marked? ☒ yes ☐ no
13. Identification number of the Notified Body 0123
14. Classification according Council Directive 93/42/EEC concerning Medical Devices Annex IX:
class: IIa rule: 2
15. Attachments, if applicable: ISO 13485 Certificate
Declaration of Conformity (DoC)
EC Certificate
List of materials
Confirmation Sterilization / Re-Sterilization
Instructions for use (IFU)

All data refer to typical single values and may be subject to alterations, valid without signature.

Edition: November 2020