

Certificate of Conformity

VWR Cleanroom Sleeve Protector

VWR International certifies that the specified products meet all VWR specifications and quality criteria for release. The products comply in every respect with the relevant particular specifications, drawings, relevant standards and regulations in force.

VWR International certifie que les produits mentionnés respectent toutes les spécifications et critères de qualité VWR. Les produits satisfont aux spécifications particulières, aux plans, aux normes et règlements en vigueur s'y rapportant.

VWR International bestätigt hiermit, daß die angegebenen Produkte alle VWR Spezifikationen und Qualitätskriterien zur Freigabe erfüllen. Die Produkte entsprechen in jeder Hinsicht den maßgeblichen Spezifikationen, den Zeichnungen sowie den relevanten Normen und gültigen Vorschriften.

VWR Catalog No: 113-7887

Lot No: 201010

Date of Manufacture : 10th Oct 2020

Date of Sterilization: 25th Oct 2020

Sterilization Expiration Date : 25th Oct 2025

Material: 60g SMS (Spunbonded Polypropylene fabric /Meltblown /Spunbonded Polypropylene fabric)

Quality System Compliance: Products are manufactured under the ISO 9001:2015 & ISO 13485:2012 standard. Products are inspected and controlled through whole production processing in accordance with current applicable product specifications and QC SOP. Inspection records are reviewed and signed off by qualified personnel for product release.

QC Testing: Inspected according to AQL level II 4.0 sampling plan per MIL-STD-105E

Break Tensile properties (N): ASTM D5034 standard, MD = 153.5, CD = 69.1

Elongation rate (%): ASTM D5034 standard, MD = 124.5, CD = 115.4

Trapezoid Tearing Strength (N): ASTM D5034 standard, MD = 159.2, CD = 82.3

Certificate of Conformity

Sterility: Sterilized by gamma irradiation, according to ISO 11137:2006, dose set meet SAL 10^{-6} with dose range of Min 27.6 kGy to Max 48.6 kGy

BSE/TSE Statement: The above product contains no ingredients of animal origin and no material derived from or exposed to animals affected by or under quarantine for Transmitting, Animal Spongiform Encephalopathy (TSE)/Bovine Spongiform Encephalopathy(BSE)

Other Characteristics: Latex Free

This document has been produced electronically and is valid without a signature

Date of Issue: **May 2021**