



ZPACK13951001 PLASTIC DRUM 5 L DIN 51 UN



SPECIFICATIONS

UN-certification :	3H1/Y1.9/200/
Material :	HDPE
Neck Type :	DIN 51
Height :	241,00 mm
Length :	182,00 mm
Wide :	162,00 mm

AVAILABLE ACCESSORIES

Key for opening :	29548.064
Tap DIN 51 :	29548.042

Remark : Packaging materials and specifications shown may be modified without preceding notice
Version dd. 22 mai 2012

tel/tél +32(0)16 385 011
fax +32(0)16 385 385
info@be.vwr.com
www.vwr.com

Researchpark Haasrode 2020
Geldenaaksebaan 464 B-3001 Leuven
RPR/RPM Leuven
BTW/TVA BE 0403.593.343