

Certificate of Quality

VWR Catalogue Number	See list
Description	VWR® Disposable Spatulas, PP
Country of Origin	Manufactured in China
Date of Issue (yyyy-mm-dd)	2026-01-07

VWR Catalogue Number	Description	Color	Length	Sterile	Unit
80081-194	Micro antistatic disposable spatula, double-ended micro tip	Opaque	14 cm (5½")	Non sterile	Pack of 300
					Case of 7500
89233-128	Eco short disposable spatula, double-ended spoon and scoop	Opaque	14 cm (5½")	Non sterile	Pack of 300
					Case of 7200
80081-188	Standard disposable spatula, double-ended spoon and scoop	Blue	21 cm (8¼")	Non sterile	Pack of 300
					Case of 7500
80081-190	Standard disposable spatula, double-ended spoon and scoop	Opaque	21 cm (8¼")	Non sterile	Pack of 300
					Case of 7500
89097-816	Sterile standard disposable spatula, double-ended spoon and scoop	Opaque	21 cm (8¼")	EO sterilisation	Pack of 100
					Case of 500
10816-246	Sterile nontoxic disposable spatula, double-ended spoon and scoop, nonpyrogenic	Opaque	21 cm (8¼")	EO sterilisation	Pack of 100
					Case of 500
80081-192	Macro disposable spatula, double-ended spoon and pick	Green	31 cm (12")	Non sterile	Pack of 150
					Case of 3600

Quality System Compliance

This certificate confirms that the product has been sourced from a supplier or manufacturer who has declared that random samples of each production batch are reviewed for adherence to specifications. SGS testing is performed annually on random samples to assure compliance with US and EU food contact regulations. Also, sterilization is done in compliance with **ISO 11135** requirements.

QC Testing

Representative products samples are collected and inspected in accordance with current applicable QC requirements.

Product Specifications

Material	Polypropylene (PP), Detailed material statements are provided in the respective batch documents of each product.
Sterilization	The products labelled sterile are sterilized using ETO (Ethylene Oxide) as the sterilant to an SAL of 10^{-6} according to the guidelines of ISO 11135.
ATP Assay	No ATP testing is required as this product is for research use only and is not a medical device.
DNase & RNase Free	The products are not tested for DNase & RNase.
BSE/TSE	No use of any raw material produced from or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE). The product 80081-194 contains ~0.4% of the AMPACET40390 Hostat10 BOPP Masterbatch which is made of one or more ingredients formulated with animal-derived substances that comply with USDA/APHIS Title 9 CFR, EMEA/410/01 (rev. 10/02/02), EU regulation 1774/2002/EC, Directive 2000/6/EC and Commission Directive 2003/32/EC.
Cytotoxicity	The products are not tested for Cytotoxicity
Non-Pyrogenic Statement	All sterile products are further processed and tested to meet Endotoxin Limits per USP 85 requirements (≤ 0.5 EU/mL or 20 EU/device). The testing is not performed on non-sterile products.
Latex Statement	The product is free from latex.
BPA Statement	Bisphenols are not used in the manufacture of raw material and is not expected to be present.
DEHP Statement	Diethylhexyl phthalate (DEHP) is not included in the manufacture or the formulation of the raw material and are not expected to be present.
RoHS	Not Applicable.
REACH Statement	No substances as listed on >Substance of Very High Concern< candidate list published are used in manufacturing of the raw materials used to produce this product. No routinely analyse is performed.

Storage Conditions	Store at Room temperature under normal lab conditions out of direct sunlight.
Shelf-Life	Products labelled as sterile have a shelf life of 72 months whereas non-sterile products have unlimited shelf life.

Disclaimer: VWR'S Quality Management System (QMS) is ISO 9001:2015 certified as a wholesale distributor of finished products, ensuring effective oversight of finished product handling, traceability and distribution processes. VWR is not a manufacturer and does not engage in manufacturing, co- manufacturing, or processing activities such as design, packaging, labelling, testing, or quality release. Consequently, the supplier or manufacturer's QMS certification may differ from VWR's.

VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. This product is intended for research use only. It is not for use in diagnostic procedures, therapeutic applications, or any human or clinical use. This document has been produced electronically and is valid without a signature.

Version History

Sl.no	Version number	Description	Month & Year
1	1.0	Created General CoQ using new layout	January 2026