



KIMTECH PREP*
KIMTEX* Wipers
Certificate of Conformance

This document certifies that KIMTECH PREP* KIMTEX* Wipers conform to Kimberly-Clark Corporation's internal specifications for this product as stated below. We use a system of in-process and audit inspections to assure that only product which meets specification is released for shipment to customers.

Certified by: *Mekelo White* Comments: SW30151

Code	Type	Color/ Colour	Material	Size	Count
33560	¼ Fold	Blue	Meltblown	12.5" x 12" 31.8 x 30.5 cm	66 wipers/pack 8 packs/case
33570	Pop-Up* Box	Blue	Meltblown	8.8" x 16.8" 22.4 x 42.7 cm	100 wipers/box 5 boxes/ case
33352	BRAG* Box	Blue	Meltblown	12.1" x 16.8" 30.7 x 42.7 cm	180 wipers/box 1 boxes/case
33241	Jumbo Roll	Blue	Meltblown	9.8" x 13.4" 24.9 x 34 cm	717 wipers/roll 1 roll/case

Physical property data is as produced, subsequent storage conditions and the age of the product, which are beyond Kimberly-Clark's control, may affect the properties of the product as received.

	Average	Test Method
CD Wet Grab		
lbs	8.5	ASTM D 5034; STM 5668
MD Grab		
lbs	8.8	STM 5668
Water Capacity		
g/4" x 4"	3.5	STM 2401
Mineral Oil Capacity		
g/4" x 4"	3.9	STM 2401

Should you require additional information, please contact our Technical Help Line at 888/346-4652.