



Certificate of Conformity

VWR Dispenser Tip 10 ml, Universal

VWR International certifies that the specified products meet all VWR specifications and quality criteria for release. The products comply in every respect with the relevant particular specifications, drawings, relevant standards and regulations in force.

VWR International certifie que les produits mentionnés respectent toutes les spécifications et critères de qualité VWR. Les produits satisfont aux spécifications particulières, aux plans, aux normes et règlements en vigueur s'y rapportant.

VWR International bestätigt hiermit, daß die angegebenen Produkte alle VWR Spezifikationen und Qualitätskriterien zur Freigabe erfüllen. Die Produkte entsprechen in jeder Hinsicht den maßgeblichen Spezifikationen, den Zeichnungen sowie den relevanten Normen und gültigen Vorschriften.

VWR Catalog No: 613-5873

Lot No: 947416/1

Date of Manufacture: June 2020

Expiry Date: June 2025

Properties: Human DNA, RNase, and ATP Free.

Sterility: Irradiated according to DIN EN ISO 11137:2006 with beta-rays with a dosage of minimum 15Kgy

All test performed by an independent certified testing laboratory.

This document has been produced electronically and is valid without a signature

Date of Issue: **January 2021**