

Product Specification Sheet

VWR BOXES FOR BROKEN GLASS

VWR International certifies that the specified products meet all VWR specifications and quality criteria for release. The products comply in every respect with the relevant particular specifications, drawings, relevant standards and regulations in force.

VWR International certifie que les produits mentionnés respectent toutes les spécifications et critères de qualité VWR. Les produits satisfont aux spécifications particulières, aux plans, aux normes et règlements en vigueur s'y rapportant.

VWR International bestätigt hiermit, daß die angegebenen Produkte alle VWR Spezifikationen und Qualitätskriterien zur Freigabe erfüllen. Die Produkte entsprechen in jeder Hinsicht den maßgeblichen Spezifikationen, den Zeichnungen sowie den relevanten Normen und gültigen Vorschriften.

Cat. No.	WxDxH(mm)	Description	Pk
119-0036	305x305x685	Floor*	6
119-0037	203x203x254	Bench**	6

* Maximum weight: 18 kg

** Maximum weight: 5.5 kg

Material: Cardboard carton with PE pouch 50 µm thick.

Application:

- Ideal for collecting, transporting and removing broken glass from the laboratory.
- Consists of a sturdy box made of strong corrugated cardboard and a PE insertion pouch which prevents any possible escape of fluids adhering to the broken glass.
- Folded flat for space saving storage
- A flap is built into the lid for easy filling
- The internal safety lid is closed for the removal of the container
- Inner Flap is pushed/pulled up to close off opening prior to disposal
- Maximum Weight of Glass Waste 12 lbs. / 5.5kg
- Does not Qualify as a DOT container

Disclaimer: VWR wants to state that these specifications can only be regarded as general information and will not discharge the user from his obligation to ensure the product is suitable for the intended use - If there is any doubt ask us.

Date of issue: **May 2013**