

Standard single use syringes

Single use syringes Luer TBC scale



1. The narrow silicone seal minimizes the contact area between the silicone and the solution so that solid particle contamination is reduced considerably. The silicone seal provides easy and accurate dosage management.
2. Luer-Slip centric nozzle. TBC scale. The dual scale provides easy and precise reading irrespective of the handling of the syringe.
3. The plunger requires despite of reduced amount of silicone oil a minimum use of force due to smooth plunger action even when used single-handed. Finger grip of an ergonomic design makes the syringe easy to handle.



Disclaimer:

The above is indicated only for information. This document is not a specification. The values given are subject to the usual fluctuations. Please follow the instructions for use by the medication manufacturer and the actual SPC (Summary of Product Characteristics). Common hygiene measures and working methods are to be followed at all times.

Subject to alterations without prior notice!

REF	Product description	Packaging	GTIN	Filling/Storage volume
62.1202	1 ML ONCE® Luer TBC 0,01 grad	1 100 3200 51200	05706225000034 05706225002618 05706225000010 -	



Standard single use syringes

Single use syringes Luer TBC scale

REF	Product description	Packaging	GTIN	Filling/Storage volume
62.1612	1 ML CODAN Luer TBC 0,01 grad	1	05706225000096	
		100	05706225001819	
		3200	05706225000072	
		51200	-	
62.1657	1 ML CODAN Luer TBC 0,01 grad	1	05706225000188	
		100	05706225001833	
		3200	05706225000164	
		51200	-	