



Certificate of Conformity

VWR®, Protection Hood and Face Cover
113-8202 / 114-9456

CHARACTERISTICS:

- White SPP material 20 g/m²
- Bulk package
- Latex free
- Size U and XL options

MATERIAL SPECIFICATION


20gsm SPP

Item		Test method	average value standard	
			average value	tolerance for average value
Weight(g/m ²)		GB/T4669-1995	20.00	+/-2
Break Tensile properties(N)	MD	ASTM D5034	>40	
	CD	ASTM D5034	>20	
Elongation rate(%)	MD	ASTM D5034	>60	
	CD	ASTM D5034	>70	



Certificate of Conformity

DIMENSION



Unit:CM

ECN	Colour	Size	Dimension								
			A	B	C	D	E	F	G	H	I
113-8202	White	Universal	34.9	38.1	21.0	3.8	14.0	14.0	17.8	24.1	15.2
114-9456	White	XL	37.5	38.1	21.0	3.8	14.0	16.5	17.8	24.1	15.2
Tolerance			±1	±1	±1	±1	±2	±1	±2	±1	±1

Remark: E/G: Stretch size

PACKAGE

Package/Quantity	100pcs/bag, 5bags/Case
Best storage Conditions	10° - 30°C
	40% - 75% humidity
	dry
Country of Origin	China



Certificate of Conformity

QUALITY ASSURANCE CRITERIA:

Lab coats are guaranteed to meet 4.0 AQL level II sampling plan per MIL-std-150E for the following defect criteria:

- | | | | |
|----|------------------------|----|------------------------------------|
| 1. | Broken elastic band | 4. | Holes/tears in fabric material |
| 2. | Irregular SPP material | 5. | Physical dimensions non-compliance |
| 3. | Dirty | 6. | Not clean |

Disclaimer: VWR wants to state that this declaration will not discharge the user from his obligation to ensure the product is suitable for the intended use. - If there is any doubt ask us.

Date of issue: **Apr 2025**