



Specification

1.07661.0500 Lactose broth acc. FDA-BAM GranuCult®

Specification

Appearance (clearness)	clear
Appearance (color)	yellowish
pH-value (25 °C)	6.7 - 7.1

Growth promotion test in accordance with the current version of DIN EN ISO 11133.

Specification

Inoculum on reference medium (Escherichia coli ATCC 8739 (WDCM 00012))	≤ 100
Inoculum on reference medium (Escherichia coli ATCC 25922 (WDCM 00013))	≤ 100
Inoculum on reference medium (Salmonella typhimurium ATCC 14028 (WDCM 00031))	≤ 100
Inoculum on reference medium (Salmonella enteritidis ATCC 13076 (WDCM 00030))	≤ 100
Inoculum on reference medium (Salmonella diarizonae ATCC 12325)	≤ 100
Inoculum on reference medium (Salmonella diarizonae ATCC 29934)	≤ 100
Inoculum on reference medium (Salmonella abortusovae ATCC 9842)	≤ 100
Inoculum on reference medium (Enterococcus faecalis ATCC 19433 (WDCM 00009))	≤ 100
Inoculum on reference medium (Enterococcus faecalis ATCC 29212 (WDCM 00087))	≤ 100
Growth (Escherichia coli ATCC 8739 (WDCM 00012))	good to very good
Growth Escherichia coli ATCC 25922 (WDCM 00013)	good to very good
Growth (Salmonella typhimurium ATCC 14028 (WDCM 00031))	good to very good
Growth (Salmonella enteritidis ATCC 13076 (WDCM 00030))	good to very good
Growth (Salmonella diarizonae ATCC 12325)	good to very good
Growth (Salmonella diarizonae ATCC 29934)	good to very good
Growth (Salmonella abortusovae ATCC 9842)	good to very good
Growth (Enterococcus faecalis ATCC 19433 (WDCM 00009))	good to very good
Growth (Enterococcus faecalis ATCC 29212 (WDCM 00087))	good to very good
Gas (Escherichia coli ATCC 8739 (WDCM 00012))	+
Gas (Escherichia coli ATCC 25922 (WDCM 00013))	+
Gas (Salmonella typhimurium ATCC 14028 (WDCM 00031))	-

Specification

1.07661.0500 Lactose broth acc. FDA-BAM GranuCult®

Gas (Salmonella enteritidis ATCC 13076 (WDCM 00030))	-
Gas (Salmonella diarizonae ATCC 12325)	-
Gas (Salmonella diarizonae ATCC 29934)	-
Gas (Salmonella abortusovine ATCC 9842)	-
Gas (Enterococcus faecalis ATCC 19433 (WDCM 00009))	-
Gas (Enterococcus faecalis ATCC 29212 (WDCM 00087))	-

Incubation: 24 ± 2 Stunden bei 35 ± 1 °C, aerob

Stefanie Fischer

Responsible laboratory manager quality control

This document has been produced electronically and is valid without a signature.