



Product Specifications sheet

VWR Collection Measuring Flasks, Clear, Class A

VWR International certifies that the specified products meet all VWR specifications and quality criteria for release. The products comply in every respect with the relevant particular specifications, drawings, relevant standards and regulations in force.

VWR International certifie que les produits mentionnés respectent toutes les spécifications et critères de qualité VWR. Les produits satisfont aux spécifications particulières, aux plans, aux normes et règlements en vigueur s'y rapportant.

VWR International bestätigt hiermit, daß die angegebenen Produkte alle VWR Spezifikationen und Qualitätskriterien zur Freigabe erfüllen. Die Produkte entsprechen in jeder Hinsicht den maßgeblichen Spezifikationen, den Zeichnungen sowie den relevanten Normen und gültigen Vorschriften.

VWR Cat. N°	Capacity (ml)	Joint Size	Tolerance (± ml)
612-3738	5	10/19	0,040
612-3740	10	10/19	0,040
612-3741	20	10/19	0,040
612-3742	25	10/19	0,040
612-3743	50	12/21	0,060
612-3744	100	14/23	0,100
612-3745	200	14/23	0,150
612-3746	250	14/23	0,150
612-3818	500	19/26	0,250
612-3819	1000	24/29	0,400
612-3820	2000	29/32	0,600

Intended uses: Laboratory applications

Material: Borosilicate 3.3 glass

Shape: Measuring Flasks Class A, clear glass

Compliant to: DIN EN ISO 4788 Borosilicate glass 3.3 - properties
ASTM E438, Type 1, Class A
Eur. Pharmacopoeia 5.0, glass type I

Chemical Properties: Have been tested by the Central Glass / Ceramic Research Institute, Kolkata, India and the report is found compliant to characteristics of Boro 3.3 Glass which should also comply to requirements of EN ISO 3585, art.4.1,4.2,4.3a, 4 and ASTM E438-92. The glass should also have following classification:
Hydrolytic class HGB 1 according to ISO 719
Hydrolytic class HGA 1 according to ISO 720



Product Specifications sheet

VWR Collection Measuring Flasks, Clear, Class A

Disclaimer: VWR wants to state that these specifications can only be regarded as general information and will not discharge the user from his obligation to ensure the product is suitable for the intended use - If there's any doubt ask us.