



Quality Certificate

VWR® STERILE BIN LINERS

As per the manufacturer, the below product meets the following criteria:

North American Catalog No:	89495-442
Lot Number:	021687-000
Description:	VWR Binliner Irradiated 40 X 48 CS50
Date of Manufacture:	September 11, 2020
Shelf Life:	36 Months

- Width: 40"
- Length: 48"

Gamma Process Run ID: 123807A

Processing Run Date: September 26, 2020

Specified Dose (kGy): Minimum 25.0; Maximum 40.0

Delivered Dose (kGy): Minimum 28.5; Maximum 38.1

Signed:

A handwritten signature in black ink that reads "Jamie Ethier".

Jamie Ethier
VP Global Quality
VWR, Part of Avantor

Date: October 16, 2020