

Certificate of Conformity

VWR Catalogue Number	732-2386
Description	PCR plates, 96-well with Coloured Lettering, Low profile, non skirted, Clear, with 0.15 mL Well Volume
Lot Number	0000005713
Country of Origin	Manufactured in Germany
Date of Issue (yyyy-mm-dd)	2022-04-05

Quality System Compliance

Products are manufactured under the **ISO 9001:2015** standard. Products are inspected and controlled through whole production processing in accordance with current applicable product specifications QC SOP and Clean room Class 8 according to ISO 14644-1.

QC Testing

Representative products samples are collected and inspected in accordance with current applicable QC requirements.

Product Specifications

Material	Polypropylene.
Sterilization	Not Applicable.
ATP Assay	Not Applicable.
Endotoxin Free	This product is free of Endotoxin.
DNase & RNase Free	This product is free of any detectable DNase/RNase contamination.
PCR Inhibitors	This product is free of PCR Inhibitors.
BSE/TSE	Tallow derived additives may be used in the manufacture of this product. "Mad Cow" tallow derived materials used in this product fulfill the requirements laid down in the Regulations 1069/2009/EC, and 142/2011/EC, and the "Note for Guidance EMA/410/01, and as amended.

Cytotoxicity	Not Applicable.
Non-Pyrogenic Statement	The acceptance level for product is 0.01 EU/mL device (LAL Gel Clot Method).
Latex Statement	This product is not made from natural rubber latex. Bisphenol are not used in the manufacture of the raw material and are not expected to be present.
DEHP Statement	Please see RoHS information below.
RoHS	No substances (Lead, Cadmium, Mercury, Hexavalent Chromium (Cr6+), Poly Brominated Biphenyls (PBB), Poly Brominated Diphenyl ethers (PBDE), Bis(2-Ethyhexyl) Phthalate (DEHP), Benzyl Butyl Phthalate (BBP), Dibutyl Phthalate (DBP), Diisobutyl Phthalate (DIBP)) are used in manufacturing the raw materials and final product. No routinely analyse is performed.
REACH Statement	No substances as listed on >Substance of Very High Concern< candidate list published 2022-01-17 are used in manufacturing of the raw materials used to produce this product. No routinely analysis is performed.
Storage Conditions	Store in a Room temperature.

Disclaimer: VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. The purpose of the product is for use in laboratory only.
This document has been produced electronically and is valid without a signature.