

Buffer Concentrate
DILUT-IT
pH 7.00 [Phosphate] Green



Material No: 4861
Effective from: 2015/06/30

Product Specification

Test	Specification
pH at 20°C	6.98 – 7.02

For Laboratory, Research or Manufacturing use.
Volumetric Concentrate, for dilution to 500 ml.
This buffer is traceable to NIST Standard Reference Buffer.
Contains Potassium Dihydrogen Phosphate / Disodiumhydrogen Phosphate / Coloring Agent.

Barbara Rykowska
Quality Assurance Coordinator

For questions on this Product Specification please contact Technical Service at technical.service@avantorsciences.com
Avantor Performance Materials Poland S.A. Sowinskiego 11 Gliwice 44-101 Poland Tel: +48 322392000