

Certificate of Conformity

VWR Catalogue Number	479-2213
Description	Freezer box, Cardboard, white with Divider
Country of Origin	Manufactured in Czech Republic
Date of Issue (yyyy-mm-dd)	2022-03-04

Quality System Compliance

Products are manufactured under the **ISO 9001:2016 & ISO 14001:2016** certified plant. Products are inspected and controlled with current applicable product specifications.

QC Testing

Representative products samples are collected and inspected in accordance with current applicable QC requirements.

Product Specifications

Material	Cardboard	ESKA 1.0 mm; White/Grey, 670 g/m2
	Cover Paper	130 g/m2 + 1/0 lamination (glass).
	Dividers	Master Karton, MBB 300 g/m2, GC2
Temperature Limitations	-140 °C.	
Sterilization	Non-Sterile	
ATP Assay	Not Applicable.	
DNase & RNase Free	Not Applicable.	
BSE/TSE	Not Applicable.	
Cytotoxicity	Not Applicable.	
Non-Pyrogenic Statement	Not Applicable.	
Latex Statement	Latex Free.	
BPA Statement	Not Applicable.	
DEHP Statement	Not Applicable.	


RoHS	Not Applicable.
REACH Statement	Not Applicable.
Storage Conditions	No Specific storage condition required.
Shelf Life	Not Applicable.

Disclaimer: VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. The purpose of the product is for use in laboratory only.
This document has been produced electronically and is valid without a signature.

Annex

Technical Drawing

479-2213



Divider 9 x 9 for 133 x 133 x 50 mm

